Profilis Paincath®

Cathéter multiperforé  
**Profilis Paincath®**

Références	P500-05	P500-10	P500-15	P500-20	P500-25	P500-30
	P500-05R	P500-10R	P500-15R	P500-20R	P500-25R	P500-30R
	P500-05T	P500-10T	P500-15T	P500-20T	P500-25T	P500-30T
	P500-05Y	P500-10Y	P500-15Y	P500-20Y	P500-25Y	P500-30Y

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Profilis Paincath®

Multiperforated catheter  
**Profilis Paincath®**

References	P500-05	P500-10	P500-15	P500-20	P500-25	P500-30
	P500-05R	P500-10R	P500-15R	P500-20R	P500-25R	P500-30R
	P500-05T	P500-10T	P500-15T	P500-20T	P500-25T	P500-30T
	P500-05Y	P500-10Y	P500-15Y	P500-20Y	P500-25Y	P500-30Y

Profilis and Profilis Paincath® are trade marks by **distriflora medical SA**

Profilis Paincath®

**Product Description :**  
The multiperforated **Profilis Paincath®** catheter is made of medical approved polyurethane material.  
OD = 1mm ; ID = 0,5mm ; Length = 50cm.  
The catheter is perforated on 05,10,15,20,25 or 30cm, distal end side closed by a titanium round end ORX.  
The other side is opened to be linked to a Tuohy Bors Adaptor.

**Indications for Use :**  
The **Profilis Paincath®** multiperforated catheter is a medical device allowing in situ post-operatively, the diffusion of analgesic drugs for the temporary treatment of pain in orthopaedic, visceral, gynaecologic and rachis surgery.

**Contraindications:**  
None, except risks in association with pain treatment.

**Precautions :**

- Carefully read and follow all instructions prior use.
- Only qualified physicians shall insert, manipulate and remove this device.
- Medical device intended for Single Patient Use. Do not Reuse.
- Caution: do not use if the package is damaged, opened or if the expiration date is passed.
- The device is supplied in a sterile package and is non-pyrogenic. Sterilized by ethylene oxide. Do not Resterilize.

**Possible complications :**

- Allergy to the analgesic drugs.
- No analgesic drugs diffusion.
- Haematomal formation.
- Air embolus oponevrotic.
- Lesion or laceration of tissues or vessels.
- Catheter break during ablation.

**Instructions for Use :**  
If you are using a protocol to rinse the catheter, follow steps 1/2/3, otherwise go directly to step 4 :

- Assemble and lock the Tuohy Bors Adaptor on the opened end side.
- Purge, through the connector, the catheter with a syringe of 10cc, filled of sterile physiological serum. Caution : do not use syringe of less than 10cc.
- Disassemble the Tuohy Bors Adaptor.
- Set up in the open wound the distal end of the multiperforated catheter, according to the selected set:
  - Standard Set with tunneler from opened wound to outside : Réf. from P500-05 to P500-30.
  - Set with Percutaneous Introducer, from open wound to outside : Réf. from P500-05R to P500-30R.
  - Set with Tuohy needle, from outside to opened wound : Réf. from P500-05T to P500-30T.
  - 2 standard Sets, allowing therapy with 2 catheters in situ linked to a Y: Réf. P500-05Y to P500-30Y.
- For the Patient comfort, cut the surplus of the exteriorized catheter, assemble and lock the Tuohy Bors Adaptor. Fix the whole set on the skin following the protocol in place.  
The unit is ready to receive the analgesic treatment through the Tuohy Bors Adaptor, connected to an elastomeric pump or any other system of diffusion (flow < at 10 ml/h).

**Retrieval :**  
After total drug diffusion, the catheter will be delicately withdrawn from the treated zone, to avoid any break. The physician will check that the catheter was entirely withdrawn, including the titanium round end ORX.  
The elimination of this device must be made in accordance to the imposed rules of safety by the medical establishment, to avoid any risk of contamination for doctors, nurses or any other people.

Manufacturer **distriflora medical SA**  
16, rue Paul Bert  
42000 SAINT-ETIENNE - France

Sales Office **distriflora medical SA**  
B.P.14 - 69970 CHAPONNAY - Lyon - France  
Tél. : +33 (0)4 72 50 43 70 - Fax : +33 (0)4 72 50 43 67  
E-mail : [distriflora@distriflora.com](mailto:distriflora@distriflora.com) - Site : [www.distriflora.com](http://www.distriflora.com)

## Instructions d'Utilisation / *Instructions for Use*

CATHETER MULTI PERFORE  
MULTI PERFORATED CATHETER  
**Profilis Paincath**<sup>®</sup>

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P500-05S	P500-05R	P500-05Y	P500-05Q
P500-10S	P500-10R	P500-10Y	P500-10Q
P500-15S	P500-15R	P500-15Y	P500-15Q
P500-20S	P500-20R	P500-20Y	P500-20Q
P500-25S	P500-25R	P500-25Y	P500-25Q
P500-30S	P500-30R	P500-30Y	P500-30Q
P1000-05S	P1000-05R	P1000-05Y	P1000-05Q
P1000-10S	P1000-10R	P1000-10Y	P1000-10Q
P1000-15S	P1000-15R	P1000-15Y	P1000-15Q
P1000-20S	P1000-20R	P1000-20Y	P1000-20Q
P1000-25S	P1000-25R	P1000-25Y	P1000-25Q
P1000-30S	P1000-30R	P1000-30Y	P1000-30Q

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CE mark obtaining date : 2008, march.

**Product Description :**

The multi perforated **PAINCATH®** catheter is made of medical approved polyurethane material. OD =1mm ; ID =0,5mm ; Length = 50 or 100 cm. The catheter is perforated of micro-holes on 05,10,15,20,25 or 30cm, distal end side closed by a titanium round end ORX. The other side is opened to be linked to a Tuohy Borst Adaptor.

**Indications for Use :**

The **Profilis Painscath®** multi perforated catheter is a medical device allowing in situ post-operatively, the diffusion of analgesic drugs for the temporary treatment of pain in orthopaedic, visceral, gynaecologic and rachis surgery.

**Contraindications:**

None, except risks in association with pain treatment.

**Precautions : CAUTION**

- Only qualified physicians shall insert, manipulate and remove this device.
- **Death warning: Neither implantation nor infusion by intravenous for analgesic.**
- **Do not Reuse.** Medical device intended for **Single Patient Use.**
- **Do not use if the package is damaged, opened or if the expiration date is passed.**
- **Do not Resterilize.** The device is supplied in a sterile package and is non-pyrogenic. Sterilized by ethylene oxide.

**Possible complications :**

- Allergy to the analgesic drugs.
- No analgesic drugs diffusion.
- Haematomas formation.
- Air embolus aponevrotic.
- Lesion or laceration of tissues or vessels
- Catheter break during ablation.

**Instructions for Use :**

If you are using a protocol to rinse the catheter, follow steps 1/2/3, otherwise go directly to step 4 :

- 1- Assemble and lock the Tuohy Borst Adaptor on the opened end side.
- 2- Purge, through the connector, the catheter with a syringe of 10cc, filled of sterile physiological serum. **Caution:** do not use syringe of less than 10cc.
- 3- Disassemble the Tuohy Borst Adaptor.
- 4- Set up in the open wound the distal end multi perforated catheter, according to the selected set:
  - a) Standard Set with tunneler from opened wound to outside : Ref. from **P500-05S** to **P500-30S** and **P1000-05S** to **P1000-30S**
  - b) Set with Percutaneous Introducer, from open wound to outside : Ref. from **P500-05R** to **P500-30R** and **P1000-05R** to **P1000-30R**
  - c) 2 standard Sets, allowing therapy with 2 catheters in situ linked to a **Y**: Ref. **P500-05Y** to **P500-30Y**, and **P1000-05Y** to **P1000-30Y**
- 5- For the Patient comfort, cut the surplus of the exteriorized catheter, assemble and lock the Tuohy Borst Adaptor. Fix the whole set on the skin following the protocol in place.

The unit is ready to receive the analgesic treatment through the Tuohy Borst Adaptor, connected to an elastomeric pump.

**Retrieval :**

After total drug diffusion, the catheter will be delicately withdrawn from the treated zone, to avoid any break. The physician will check that the catheter was entirely withdrawn (distal end marked).

The elimination of this device must be made in accordance to the imposed rules of safety by the medical establishment, to avoid any risk of contamination for doctors, nurses or any other people.

References ending by **Q** are for box containing 10 sets, including for each set :

- 1- multi perforated **Painscath®** cathéter
- 1- Tuohy Borst Adaptor Luer lock.

Fabricant / Manufacturer :

districlass médical SA

16 Rue Paul Bert,  
42000 SAINT-ETIENNE - France

Bureaux commerciaux / Sales Office :

districlass médical SA

B.P. 14, 69970 CHAPONNAY – France  
Tel. +33 (0) 4 72 50 43 75  
Fax. +33 (0) 4 72 50 43 67  
Website : [www.districtclass.com](http://www.districtclass.com)  
e-mail : [districtclass@districtclass.com](mailto:districtclass@districtclass.com)

Fabricant / Manufacturer :

districlass médical SA

16 Rue Paul Bert,  
42000 SAINT-ETIENNE - France

Bureaux commerciaux / Sales Office :

districlass médical SA

B.P. 14, 69970 CHAPONNAY – France

Tel. +33 (0) 4 72 50 43 75

Fax. +33 (0) 4 72 50 43 67

e-mail : [districlass@districlass.com](mailto:districlass@districlass.com)

Web : [www.districtclass.com](http://www.districtclass.com)

Noti004-02 Rev2010-06

districlass médical SA

## Instructions d'Utilisation / *Instructions for Use*

### CATHETER MULTIPERFORE MULTIPERFORATED CATHETER **Profilis Paincath®**

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districtclass médical SA

**Carefully, read and follow all instructions prior use.**

Medical device **class IIb** for short term use (< 29 days).  
Warning : The maximum running time with our 0.2µm filter is 96 hours.

**Product Description:**

The multiperforated **Profilis Paincath**® catheter is made of medical approved polyurethane material.  
OD =1mm; ID =0,5mm;;  
Length = 50cm or 100cm

The catheter is perforated with micro holes on 02.5, 05, 07.5, 10, 15, 20, 25 or 30cm, distal end side closed.

The other catheter side is opened to allow the connection on its already assembled connector or on its removable Tuohy Borst connector.

**Indications for use:**

The **Profilis Paincath**® multiperforated catheter is a medical device allowing the diffusion in situ post-operatively, of analgesic drugs for the temporary treatment of pain in orthopaedic, visceral, gynaecologic and rachis surgery.

**Contraindications:**

None, except risks in association with the pain treatment.

**Precautions : WARNING**

- Only qualified physicians shall insert, manipulate and remove this device.
- **Death warning: Neither implantation nor diffusion by intravenous for analgesics.**
- **Do not reuse. Medical device intended for Single Patient Use.**
- **Do not use if the package is damaged, opened or if the expiry date is passed.**
- **Do not resterilize.**  
This device is supplied in a sterile package and is non-pyrogenic. Sterilized by ethylene oxide.  
On no account these products have to be resterilized and reused, entailing a biological major risk for the patient such as bacteria, virus, infection...

**Possible complications :**

- Allergy to the analgesic drugs.
- No analgesic drug diffusion.
- Haematomas formation.
- Aponevrotic air embolus.
- Lesion or laceration of tissues or vessels.
- Catheter breaks during ablation.

**Instructions for use:**

1- Set up in opened wound **ALL** the catheter perforation length (from the closed distal end to the marking position of the last perforation), according to the selected set:

- a) Technique with tunneler, from opened wound to the outside: Ref. **S, QS.**
- b) Technique with peel-away introducer needle, from the outside to the opened wound: Ref. **R, QR.**
- c) Technique for the treatment with 2 catheters in situ linked to a **Y**, with tunneler or peel-away introducer needle: Ref. **Y.** (accessory provided on request).

2- Assemble and lock the removable Tuohy Borst connector for the sets: Ref. **S, QS.**

3- Check the perfect catheter functioning by using a 10cc syringe filled with sterile physiological serum.

(**Caution**, preferably, use a 10cc syringe).

4- Close the wound.

5- Fix the whole set on the skin to prevent any catheter movement, following the protocol in place.

(**Caution**, do not use suture to prevent kinking of the catheter. Prefer the use of a skin adhesive Patch).


6- Connect the catheter to its elastomeric pump, or any other diffusion system (Flow < 10 ml/h) to administrate the analgesic treatment. Un filter can be dispatched upon request.

**Retrieval:**

After total drug diffusion, the catheter will be delicately withdrawn from the treated zone, to avoid any break. The physician will check that the catheter is entirely withdrawn (distal end marked).

This device disposal must be made in accordance to the imposed rules of safety by the medical establishment, to avoid any risk of contamination for doctors, nurses or any other people.

References starting by letter **Q** are for boxes containing 10 sets (See the set composition on labelling).

districlass médical SA

Instructions d'Utilisation / *Instructions for Use*

**CATHETER MULTIPERFORE**  
*MULTIPERFORATED CATHETER*

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districlass médical SA


**Carefully, read and follow all instructions prior use.**

**Instructions for use:**


Fabricant / Manufacturer :

districlass médical SA

16 Rue Paul Bert,  
42000 SAINT-ETIENNE -  
France

Bureaux commerciaux / Sales Office :

districlass médical SA

B.P. 14, 69970 CHAPONNAY – France


Tel. +33 (0) 4 72 50 43 75

Fax. +33 (0) 4 72 50 43 67

e-mail : [districlass@districlass.com](mailto:districlass@districlass.com)

Web : [www.districtclass.com](http://www.districtclass.com)

districlass médical SA

Instructions d'Utilisation / *Instructions for Use*

**CATHETER MULTIPERFORE**  
**MULTIPERFORATED CATHETER**

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**Read and follow carefully all instructions before use.**

**Product description:**

The multiperforated **Profilis Paincath** catheter is made of medical approved polyurethane material.

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**Indications for use:**

The multiperforated catheter **Profilis Paincath** is a medical device allowing the diffusion in situ post-operatively, of analgesic drugs for the temporary treatment of pain in orthopaedic, visceral, gynaecologic and rachis surgery.

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**Precautions: WARNING**

- Only qualified physicians shall insert, manipulate and remove this device.
- **Death warning: Neither implantation nor diffusion by intravenous for analgesics.**
- **Do not reuse. Single use medical device.**
- **Do not use if the package is damaged, opened, or if the date is passed.**
- **Do not re-sterilize.**

This device is supplied in a sterile package.

Sterilized by ethylene oxide.

In any case these products shall be resterilized and reused, entailing a biological major risk for the patient such as bacteria, virus, infection...

**Possible complications:**

- Allergy to the analgesic drugs.
- No analgesic drug diffusion.
- Haematomas formation.
- Aponevrotic air embolus.
- Lesion or laceration of tissus or vessels.
- Catheter breaks during ablation.

**Instructions for use:**

- 1- Set up in opened wound **ALL** the catheter perforation length (from the closed distal end to the marking position of the last perforation), according to the selected set:
  - a) Technique with tunneler, from opened wound to the outside: Ref. **S, QS, QP** (shorter for pediatrics)
  - b) Technique with peel-away introducer needle, from the outside to the opened wound: Ref. **R, QR.**
  - c) Technique for the treatment with 2 catheters in situ connected by **Y**, with tunneler or peel-away introducer needle  
(Accessories provided upon request).
- 2- Assemble and lock the removable Tuohy Borst connector for the sets: Ref. **S, QS.**
- 3- Check the perfect catheter functioning by using a 10cc syringe filled with sterile physiological serum.

**Caution**, preferably, use a 10cc syringe.

- 4- Close the wound.
- 5- Fix the whole set on the skin to prevent any catheter movement, following the protocol in place.

**Caution**, do not use suture to prevent kinking of the catheter. Prefer the use of a skin adhesive Patch.

- 6- Connect the catheter to its elastomeric pump, or any other diffusion system (Flow < 10 ml/h) to administrate the analgesic treatment. A filter can be provided upon request.

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
The discard of the device and its accessories must follow the protocol of the discard of medical wastes **DASRI (Infections Risks Care Activities Waste)** related to the centers, and must not endanger the safety or health of the patient, the users or any other person until complete destruction.

References starting by letter **Q**, are for boxes containing 10 sets (See the set composition on labelling).

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

### • Product description

The multiperforated **Profilis Poinçath** catheter is made of medical approved polyurethane material.  
OD = 1mm; ID = 0,5mm;  
Length = 50cm or 100cm.

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### • Indications for use

The multiperforated **Profilis Poinçath** catheter is a medical device allowing the diffusion in situ post-operatively, of analgesic drugs for the temporary treatment of pain.

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
**YES**  
Push and turn the catheter into the tunneler

**NO**  
Do not push the tunneler on the catheter

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 **Read and follow carefully all instructions before use.**

#### • Product description

The multiperforated **Profilis PaincAth**® catheter is made of medical approved polyurethane material.  
OD = 1mm; ID = 0,5mm;  
Length = 50cm or 100cm.

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#### • Indications for use

The multiperforated **Profilis PaincAth**® catheter is a medical device allowing the diffusion in situ post-operatively, of analgesic drugs for the temporary treatment of pain.

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#### WARNING / Precautions

 **Only qualified physicians shall insert, manipulate and remove this device.**

- **Death warning: Neither implantation nor diffusion by intravenous for analgesics.**
- **Do not reuse. Single use medical device.**
- **Do not use if the package is damaged, opened after the expiry date is passed.**
- **Do not re-sterilize.**

This device is supplied in a sterile package. Sterilized by ethylene oxide. In any case these products have to be resterilized and reused, entailing a biological major risk for the patient such as bacteria, virus, infection...

#### • Possible complications


- Allergy to the analgesic drugs
- No analgesic drug diffusion
- Haematomas formation
- Aponevrotic air embolus
- Lesion or laceration of tissue or vessels
- Catheter breaks during ablation

#### • Instructions for use

1 - Set up in opened wound **ALL** the catheter perforation length (from the closed distal end to the marking position of the last perforation), according to the selected set:

a) Technique with tunneler: from opened wound to the outside / Ref. **S, QS, QP** (shorter for pediatrics).

 **Connection of the catheter on the tunneler**

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**NO**

Do not push the catheter on the tunneler

b) Technique with peel-away introducer needle: from the outside to the opened wound / Ref. **R, QR**.

c) Technique for the treatment with 2 catheters in situ linked to a **Y**, with tunneler or peel-away introducer needle (Accessories provided upon request).

2 - Assemble and lock the removable Tuohy Borst connector for the sets / Ref. **S, QS**.

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 **Caution, preferably, use a 10cc syringe.**

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 **Caution, do not use suture to prevent kinking of the catheter. Prefer the use of a skin adhesive patch.**

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The discard of the device and its accessories must follow the protocol of the discard of medical wastes **DASRI** (Infections Risks Care Activities Waste) related to the centers, and must not endanger the safety or health of the patient, the users or any other person until complete destruction.

References starting by letter **Q** do mean that boxes are containing 10 sets (See the set composition on labelling).

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
**USER MANUAL**

**VASCULAR ACCESS CATHETER IMPLANTABLE PORT**

 *read carefully before implantation*

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## PRODUCT DESCRIPTION

The multiperforated Profilis PaincAth® catheter is made of medical approved polyurethane material.

OD = 1mm; ID = 0,5mm;

Length = 50cm or 100cm.

The catheter is perforated with micro holes on 02.5,

05, 07.5, 10, 15, 20, 25 or 30cm, distal end side closed.

The other catheter extremity is opened to allow the connection on its pre-assembled connector or on its removable Tuohy Borst connector.

## COMPOSITION OF THE KIT

A Profilis PaincAth® kit is supplied sterile and disposable.

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### 1) INDICATIONS FOR USE

The multiperforated Profilis PaincAth® catheter is a medical device allowing the diffusion in situ post-operatively, of analgesic drugs for the temporary treatment of pain.

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### 3) POSSIBLE COMPLICATIONS

- Allergy to the analgesic drugs
- No analgesic drug diffusion
- Hematomas formation
- Aponevrotic air embolus
- Lesion or laceration of tissues or vessels
- Catheter breaks during ablation

## STERILISATION

The Profilis PaincAth® kit is delivered in sterile packaging.

 Sterilized with ethylene oxide.

The reuse of single-use devices exposes the patient or the user to risks of infection. Contamination of the device can cause injury, illness or even death of the patient. Do not re-sterilize the catheter regardless of the procedure.

 Don't re-sterilize.

Reprocessing may compromise the integrity of the device or cause failure of the device.

**The manufacturer cannot be held liable for damage caused by the reuse of the medical device or its accessories.**

The contents of the package are sterile, as long as it is neither opened nor damaged.


 Don't use if the packaging is damaged.

 Don't use this material after the expiry date stated on the packaging.

## INSTRUCTIONS FOR USE

1 - Set up in opened wound ALL the catheter perforation length (from the closed distal end to the marking position of the last perforation), according to the selected set:

a) Technique with tunneler: from opened wound to the outside / Ref. S, QS, QP (tunneler shorter for pediatrics).

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**NO**

Do no push the catheter on the tunneler

b) Technique with peel-away introducer needle: from the outside to the opened wound / Ref. R, QR.

c) Technique for the treatment with 2 catheters in situ linked to a Y, with tunneler or peel-away introducer needle (Accessories provided upon request).

2 - Assemble and lock the removable Tuohy Borst connector for the sets / Ref. S, QS, QP.

3 - Check the perfect catheter functioning by using a 10cc syringe filled with sterile physiological serum.

 you must use a 10cc syringe

4 - Close the wound.

5 - Fix the whole set on the skin to prevent any catheter movement, following the protocol in place.

 do not use suture to prevent kinking of the catheter. Prefer the use of a skin adhesive patch.

6 - Connect the catheter to its elastomeric pump, or any other diffusion system following the protocol in place.

## RETRIEVAL

After total drug diffusion, the catheter will be delicately withdrawn from the treated area, to avoid any break. The physician will check that the catheter is entirely withdrawn (distal end marked).

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References starting by letter Q do mean that boxes are containing 10 sets (See the set composition on labelling).

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## INTRODUCTION

The multiperforated **Profilis PaincAth®** catheter is a medical device allowing the diffusion in situ of analgesic drugs for post-op. pain therapy.

## DESCRIPTION AND PERFORMANCES OF THE DEVICE

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Diffusion area

- Connection to an analgesic infusion line using its Luer Lock connector.
- Diffusion of the analgesia in situ up to 48 hours.

- This kind of analgesia is an alternative to the use of morphine (side effects).
- Improvement of patient care post-op pain.
- X-Ray and NMR compatible.

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## COMPOSITION OF THE SETS

The Profilis PaincAth® set is supplied sterile, single use.

Do not re-use

## COMPOSITION OF THE S OR R SET:

- 1 X-Ray multiperforated catheter: ID  $\varnothing = 1 \text{ mm}$  ; OD  $\varnothing = 0,5 \text{ mm}$  ; Length : S version (L = 50 cm) or P version (L=100 cm)

- 1 Tunneler : S version (L= 210mm) or P version (L=100mm)

- 1 removable Luer Lock Connector

- 1 Luer Lock Syringe 10cc

## COMPOSITION OF THE R SET:

- 1 X-Ray multiperforated catheter with its Luer Lock connector: ID  $\varnothing = 1 \text{ mm}$  ; OD  $\varnothing = 0,5 \text{ mm}$  ; L = 50 cm or 100 cm

- 1 Peel-away introducer needle: 18 Gauge; Length = 10 cm

Each set R, S or P is supplied with a sterile securement device: Ref. A004P 16/17G LOCKIT PLUS™, under its own EC mark.

Upon request, in its separated sterile packaging: 1 Filter 0.2 $\mu\text{m}$  Ref. A004F

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Possibility of a simultaneous diffusion on 2 catheters with the use of a Y Luer connector – Accessory not supplied.

## INDICATIONS AND CONTRAINDICATIONS

### 1) INDICATIONS

Post-op. pain therapy

### 2) CONTRE-INDICATIONS

- Associated risks to the pain therapy
- Proximity infection
- Children less than 10 Kg.

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Do not re-sterilize

The re-use will compromise the device integrity or will lead to the device default.

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Do not use if conditioning is damaged

Do not use this device after the expiry date mentioned on the conditioning.

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Never set up this device by the intravenous way for the diffusion of analgesia. **Threatening risk.**


1 - Check the perfect catheter functioning by using a 10cc syringe, filled with sterile physiological serum.

**NEVER USE A SYRINGE INFERIOR TO 10CC.**  
Risk of burst of the catheter by overpressure.

a) **Technique with tunneler:** from the surgical cavity to the outside / S, QS ; P, QP versions

Recommendations for the connection of the catheter on the tunneler

**YES**

Push the tunneler into the catheter by a screwing movement

**NO**

**Do not push the catheter on the tunneler**

b) **Technique with peel-away introducer needle:** from the outside to the surgical cavity / R, QR versions.

2 - Set up in the opened surgical cavity ALL the perforation length (closed distal end up to the marking: position of the last perforation), according to the selected set.

3 - Close the surgical cavity.

4- Fix the catheter using the fixing device provided.

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5 - Connect the catheter to the perfusion line.

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S = Technique with tunneler

R = Technique with peel-away introducer needle

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
EN

**INSTRUCTIONS FOR USE**  
MULTIPERFORATED CATHETER

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# READ CAREFULLY BEFORE USE

## 1 - INTRODUCTION

The multiperforated **Proflis PaicAth®** catheter is a medical device allowing the diffusion in situ of analgesic drugs for post-op. pain therapy.

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The multiperforated **Proflis PaicAth®** catheter is made of biocompatible rigid opaque polyurethane.

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Length = 50 cm or 100 cm.

The catheter is perforated with micro holes on 02.5 / 05 / 07.5 / 10 / 15 / 20 / 25 or 30cm, distal end side closed.

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
## 3 - COMPOSITION OF THE SETS

The **Proflis PaicAth®** set is supplied sterile, single use.

Do not re-use.

### COMPOSITION OF THE S OR P SET

• 1 X-Ray multiperforated catheter: ID  $\varnothing = 1 \text{ mm}$  ; OD  $\varnothing = 0,5 \text{ mm}$  ; Length : L = 50 cm or 100 cm.

• 1 Tunneler : S version (L= 210mm) or P version (L=100mm)

• 1 removable Luer Lock Connector

• 1 Luer Lock Syringe 10cc

### COMPOSITION OF THE R SET

• 1 X-Ray multiperforated catheter with its Luer Lock connector : ID  $\varnothing = 1 \text{ mm}$  ; OD  $\varnothing = 0,5 \text{ mm}$  ; L = 50 cm or 100 cm

• 1 Peel-away introducer needle

Ref. AIPM4  
18 Gauge; Length = 10 cm

• 1 sterile securement device, supplied under its own EC mark in each set

Ref. A004P 16/17G LOCKIT PLUS™

### UPON REQUEST

• 1 filter 0.2 $\mu\text{m}$ . Ref. A004F.  
In its separated sterile packaging

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Post-op. pain therapy.

### 2) Contre-indications

- Associated risks to the pain therapy.
- Proximity infection.
- Children less than 10 Kg.

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## 5 - STERILIZATION

The **Proflis PaicAth®** set is supplied sterile in double packaging.

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Expiry date: 4 years and 11 months.

## 6 - INSTRUCTIONS FOR IMPLANTATION


**⚠ Only experimented people can set up this device.**

**⚠ Never set up this device by the intravenous way for the diffusion of analgesia. Death risk.**

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2. Pull the tunneler with its catheter outside the surgical cavity up to the marking of the catheter: position of the last perforation.

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**⚠ Do not screw excessively the Luer connections. Fissure risk.**

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**⚠ Never use a syringe inferior to 10 CC. Risk of burst of the catheter by overpressure.**

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1. From the outside to the surgical cavity.

2. Set up in the opened surgical cavity **ALL** the perforation length (closed distal end up to the marking of the catheter: position of the last perforation).

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**⚠ Never use a syringe inferior to 10 CC. Risk of burst of the catheter by overpressure.**

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\* Close the surgical cavity.

\* Fix the catheter using the sterile fixing device provided.

**⚠ Do not use suture for fixing the catheter.**

\* Connect the catheter to the perfusion line.

**⚠ Do not screw excessively the Luer connections. Fissure risk.**

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**⚠ The physician will make sure that the catheter has been entirely removed (distal end marked).**

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
EN

INSTRUCTIONS FOR USE  
MULTIPERFORATED CATHETER

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# READ CAREFULLY BEFORE USE

## 1 - INTRODUCTION

The multiperforated **Profilis PaincAth**® catheter is a medical device allowing the diffusion in situ of analgesic drugs for post-op. pain therapy. For the following molecules:

- 1- Bupivacaïne (Bocié et AL) , ( clinical studies : 2 ;7 ;11)
- 2- Ropivacaïne ( clinical studies : 8 ;10 ;13)
- 3- Lévobupivacaïne ( clinical studies : 12)

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The multiperforated **Profilis PaincAth**® catheter is made of biocompatible radio opaque polyurethane.

**Dimensions :** ID. Ø = 1 mm ; OD Ø = 0,5 mm ; Length = 50 cm or 100 cm.

The catheter is perforated with micro holes on 02.5 / 05 / 07.5 / 10 / 15 / 20 / 25 or 30 cm, distal end side closed.

The proximal extremity allows the connection to the perfusion line either by a pre-assembled connector or by a removable connector.

Diffusion area

- Connection to an analgesic infusion line using its Luer Lock connector.
- Diffusion of the analgesia in situ up to 48 hours.
- This kind of analgesia is an alternative to the use of morphine (side effects).
- Improvement of patient car post-op pain.
- X-Ray and NMR compatible.

## 2 techniques of implantation :

\*References with S, QS or P, QP : tunneler technique

\*References with R, QR : peel-away introducer needle technique

## 3 - COMPOSITION OF THE SETS

The **Profilis PaincAth**® set is supplied sterile, single use.

Do not re-use.

## COMPOSITION OF THE S OR P SET

• 1 X-Ray multiperforated catheter: ID Ø = 1 mm ; OD Ø = 0,5 mm ; Length : L = 50 cm or 100 cm.

• 1 Tunneler : S version (L= 210mm) or P version (L=100mm)

• 1 removable Luer Lock Connector

• 1 Luer Lock Syringe 10cc

## COMPOSITION OF THE R SET

• 1 X-Ray multiperforated catheter with its Luer Lock connector : ID Ø = 1 mm ; OD Ø = 0,5 mm ; L = 50 cm or 100 cm

• 1 Peel-away introducer needle  
Réf. AIPM4  
18 Gauge; Length = 10 cm

• 1 sterile securement device, supplied under its own EC mark in each set  
Ref. A004P 16/17G LOCKIT PLUS™

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## 4 - INDICATIONS AND CONTRAINDICATIONS

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### 2) Contre-indications

- Associated risks to the pain therapy.
- Proximity infection.
- Children less than 10Kg.

## 3) Possible complications

- Hematoma
- Air embolism aponevrotic
- Lesion or laceration of tissues or vessels
- Catheter break in situ during ablation

## 5 - STERILIZATION

The **Profilis PaincAth**® set is supplied sterile in double packaging.

Sterilization by Ethylene Oxide.

The re-use of the single use device exposes the patient or the user to infection risks. The device contamination can lead to lesions, disease, up to the patient death. Do not re-sterilize the catheter whatever the process.

Do not re-sterilize

The re-use will compromise the device integrity or will lead to the device default. The manufacturer cannot be liable of any damage coming from the re-use of the medical device or any of its accessories.

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
Never set up this device by the intravenous way for the diffusion of analgesia.  
**Death risk.**

2 techniques of implantation for the **Profilis PaincAth**® catheter :

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2. Pull the tunneler with its catheter outside the surgical cavity up to the marking of the catheter: position of the last perforation.
3. Connect and secure the removable Luer connector on the catheter.

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Never use a syringe inferior to 10 CC. Risk of burst of the catheter by overpressure.

### b) Technique with peel-away introducer needle (R, QR versions) :

1. From the outside to the surgical cavity.
2. Set up in the opened surgical cavity ALL the perforation length (closed distal end up to the marking of the catheter: position of the last perforation).
3. Rinse the catheter with a sterile physiological solution using a syringe.

Never use a syringe inferior to 10 CC.  
Risk of burst of the catheter by overpressure.

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- \* Close the surgical cavity.
- \* Fix the catheter using the sterile fixing device provided.

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The physician will make sure that the catheter has been entirely removed (distal end marked).

The discard of the device and its accessories must follow the protocol of the discard of medical wastes DASRI (Infections Risks Care Activities Waste) related to the centers, and must not endanger the safety or health of the patient, the users or any other person until complete destruction.

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
EN

**INSTRUCTIONS FOR USE**  
MULTIPEFORATED CATHETER

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## 2 - DESCRIPTION AND PERFORMANCES OF THE DEVICE

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The catheter is perforated with micro holes on 02.5 / 05 / 07.5 / 10 / 15 / 20 / 25 or 30 cm, distal end side closed up to the 3 dashes.

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### Diffusion area

- The proximal extremity allows the connection to the perfusion line either by a pre-assembled connector or by a removable connector.
- Connection to an analgesic infusion line using its Luer Lock connector.
- Diffusion of the analgesia in situ up to 48 hours maximum.
- This kind of analgesia is an alternative to the use of morphine (side effects).
- Improvement of patient care post-op pain (Scientific references: 6 ; 7 ; 8).

- X-Ray and NMR compatible 

 Do not use the tunneler or introducer needle during MR imaging.

- Free latex  Free phtalate 

### 2 techniques of implantation :

- \* References with S, Q5 or P, QP : tunneler technique
- \* References with R, QR : peel-away introducer needle technique

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-  Do not re-use.

 The re-use of the single use device exposes the patient to infection risks.

## COMPOSITION OF THE S OR P SET

- 1 X-Ray multiperforated catheter: OD = 1 mm ; ID = 0.5 mm ; Length : L = 50 cm or 100 cm.

- 1 Tunneler : S version (L = 210 mm) or P version (L = 100mm)

- 1 Removable Luer Lock connector

- 1 Luer Lock syringe of 10cc

## COMPOSITION OF THE R SET

- 1 X-Ray multiperforated catheter: OD  $\varnothing = 1\text{ mm}$  ; ID  $\varnothing = 0.5\text{ mm}$  ; Length : L = 50 cm or 100 cm.

- 1 Peel-away introducer needle

Ref. A1PM4  
18 Gauges ; Length = 10 cm

- 1 sterile securement device,

supplied under its own EC mark in each set. Ref. A004P 16/17G LOCKIT PLUS™

The references starting with Q (Q5/QP/QR) mean that the box is containing 10 sets (see components of each set on the labels).

Possibility of a simultaneous diffusion on 2 catheters with the use of a Y Luer connector - Accessory not supplied.

For safe administration, a 0,2 $\mu\text{m}$  antibacterial filter should be used on the infusion line. Accessory not supplied.

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### 1) Indications

Post-op. pain therapy for a maximum of 48 h.

### 2) Contre-indications

- Lack of consent by the patient.
- Associated risks to the pain therapy.
- Proximity infection.
- Children less than 10 Kg.
- Use in epidural anesthesia.
- Known neurological abnormalities.

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 Do not use if conditioning is damaged.

 Infectious risk.

 Expiry date: 4 years and 11 months.

## 6 - INSTRUCTIONS FOR IMPLANTATION

 **Only experimented people can set up this device.**


 Never set up this device by the intravenous way for the diffusion of analgesia. **Death risk.**

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 **Recommendations**

**YES**  
Push the tunneler into the catheter by a screwing movement

**NO**  
Do not push the catheter on the tunneler

2. Pull the tunneler with its catheter outside the surgical cavity up to the marking of the catheter: position of the last perforation

3. Connect and secure the removable Luer connector on the catheter.

 Do not screw excessively the Luer connections. Fissure risk.

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b) **Technique with peel-away introducer needle (R, QR versions) :**

1. From the outside of the surgical cavity.

2. Set up in the opened surgical cavity ALL the perforation length (closed distal end up to the marking of the catheter: position of the last perforation).

3. Rinse the catheter with a sterile physiological solution using a syringe.

For the 2 techniques of implantation :

 For rinsing, never use a syringe inferior to 10 CC.

Risk of burst of the catheter by overpressure.

- \* Close the surgical cavity.

- \* Fix the catheter using the sterile fixing device provided.

 Do not use suture for fixing the catheter.

Risk of catheter pinch-off.

- \* Connect the catheter to the perfusion line.

 Always take care of the asepsis of the injection site.

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## 7 - ELIMINATION OF DEVICE AND ACCESSORIES

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 In case of inefficiency of pain treatment with our **Profilis PaincAth®** catheter and the analgesic, a substitution treatment chosen by the medical profession will have to be performed on the patient. The catheter will be removed.

 Any serious incident involving the device will be notified to the manufacturer and the competent authority of the member state in which the user and / or patient is established.

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**1 - INTRODUCTION**

The multiperforated **Profilis PaincAth®** catheter is a medical device allowing the diffusion in situ of analgesic drugs for post-op. pain therapy. For the following molecules:

1. Bupivacaine (scientific references 1 ; 2 ; 3 ; 4)
2. Ropivacaine (scientific references 1 ; 2 ; 3 ; 4 ; 5 ; 6)
3. Lévocabupivacaine (scientific references 1 ; 2 ; 5)

**2 - DESCRIPTION AND PERFORMANCES OF THE DEVICE**

The multiperforated **Profilis PaincAth®** catheter is made of biocompatible radio opaque polyurethane.

Outside Ø = 1 mm ; Inside Ø = 0,5 mm ; Length = 50 cm or 100 cm.  
The catheter is perforated with micro holes on 02.5 / 05 / 07.5 / 10 / 15 / 20 / 25 or 30cm, distal end side closed up to the 3 dashes  
(See commercial references on a table at the end of the IFU).

**Diffusion area**

- The proximal extremity allows the connection to the perfusion line either by a pre-assembled connector or by a removable connector.
- Connection to an analgesic infusion line using its Luer Lock connector.
- Diffusion of the analgesia in situ **up to 48 hours maximum**

Biological safety risk:  
Possible local and/or toxic effects on patient.

- This kind of analgesia is an alternative to the use of morphine (side effects).
- Improvement of patient care post-op pain (Scientific references: 6; 7; 8).
- X-Ray and MRI compatible. **MR**

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- 1 Peel-away introducer needle Ref. AIPM4 18 Gauges ; Length = 10 cm

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**1) Indications**

Post-op. pain therapy for a maximum of 48 h (See risk §2).

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**3) Possible complications**

- Hematoma or allergic reactions.
- Air embolism aponevrotic.
- Lesion or laceration of tissues or vessels.
- Catheter pinch-off and reduced or absence of flow.
- Catheter break in situ during ablation.

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**STERILEEO** Sterilization by Ethylene Oxide.

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
2 implantation techniques for the catheter **Profilis PaincAth®** :

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- a) Close the surgical cavity.
- b) Fix the catheter using the sterile fixing device provided.

Do not use suture for fixing the catheter. Risk of catheter pinch-off.

- c) Connect the catheter to the perfusion line.

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In case of inefficiency of pain treatment with our Profilis PaincAth catheter and the analgesic, a substitution treatment chosen by the medical profession will have to be performed on the patient. The catheter will be removed

Any serious incident involving the device will be notified to the manufacturer and the competent authority of the member state in which the user and / or patient is established;

Kit unitaire / Unitary set		Kit par 10 / Box of 10 sets	
P500-02.5 S/R/P		P500-02.5 QS/QR/QP	
P500-05 S/R/P	P1000-05 S/R	P500-05 QS/QR/QP	P1000-05 QS/QR
P500-07.5 S/R/P		P500-07.5 QS/QR/QP	
P500-10 S/R	P1000-10 S/R	P500-10 QS/QR	P1000-10 QS/QR
P500-15 S/R	P1000-15 S/R	P500-15 QS/QR	P1000-15 QS/QR
P500-20 S/R	P1000-20 S/R	P500-20 QS/QR	P1000-20 QS/QR
P500-25 S/R	P1000-25 S/R	P500-25 QS/QR	P1000-25 QS/QR
P500-30 S/R	P1000-30 S/R	P500-30 QS/QR	P1000-30 QS/QR

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**1 - INTRODUCTION**

The multiperforated **Profilis PaincAth®** catheter is a medical device allowing the diffusion in situ of analgesic drugs for the treatment of patient pain in post-op, dedicated to surgeries such as: rachis, pediatric, orthopedic, laparotomy, urology, visceral, gynecology, caesarian, ....

The multiperforated **Profilis PaincAth®** catheter is made of biocompatible radio opaque polyurethane, free latex and free phthalate

Outside Ø = 1 mm; Inside Ø = 0,5 mm; Length = 50 cm or 100 cm.  
The catheter is perforated with micro holes on 2.5 / 5 / 7.5 / 10 / 15 / 20 / 25 or 30cm, distal end side closed up to the 3 dashes.

(See commercial references on a table at the end of the IFU).

TARGET POPULATION: Patient with a weight over 10 kg without age limit.

**2 - PRODUCT PERFORMANCES AND BENEFICES**

**>Clinical Performances :**

- The multiperforated catheter **Profilis PaincAth®** is a medical device allowing the in situ diffusion of analgesic drugs up to **48h00 maximum** for the treatment of post op pain for the following molecules:

1. Bupivacaine (scientific references 1 ; 2 ; 3 ; 4)
2. Ropivacaine (scientific references 1 ; 2 ; 3 ; 4 ; 5 ; 6)
3. Lévocabupivacaine (scientific references 1 ; 2 ; 5)

**>Clinical Benefices :**

- Continuous Analgesia effective for 48 hours (Scientific References: 6 ; 7 ; 8).
- Alternative to opioid consumption and their associated side effects
- Faster post-operative recovery;
- Decrease in the hospitalization length.

**>Mechanical Performances :**

- The proximal extremity allows the connection to the perfusion line either by a pre-assembled connector or by a removable connector.
- Connection to an analgesic infusion line using its Luer Lock connector.

• X-Ray and MRI compatible.

Do not use the tunneler or introducer needle during MR imaging.

**3 - COMPOSITION OF THE SETS**

The **Profilis PaincAth®** set is supplied sterile, single use.

Do not re-use.

The re-use of the single use device exposes the patient to infection risks.

**2 techniques of implantation:**

- \* References with S, QS or P, QP: tunneler technique
- \* References with R, QR: peel-away introducer needle technique

**COMPOSITION OF THE S OR P SET**

- 1 X-Ray multiperforated catheter: OD = 1 mm ; ID = 0,5 mm ; Length : L = 50 cm or 100 cm.
- 1 Tunneler : S version (L = 210 mm) or P version (L = 100mm)

- 1 Removable Luer Lock connector

- 1 Luer Lock syringe of 10cc

**COMPOSITION OF THE R SET**

- 1 X-Ray multiperforated catheter: OD Ø = 1 mm ; ID Ø = 0,5 mm ; Length : L = 50 cm or 100 cm.

- 1 Peel-away introducer needle  
Ref. AIPM4  
18 Gauges ; Length = 10 cm

- 1 sterile securement device, supplied under its own EC mark in each set. Ref. A004P 16/17G LOCKIT PLUS™

The references starting with Q (QS/QP/QR) mean that the box is containing 10 sets (see components of each set on the labels).

Possibility of a simultaneous diffusion on 2 catheters with the use of a Y Luer connector. Accessory not supplied.

For safe administration, a 0,2µm antibacterial filter should be used on the infusion line. Accessory not supplied.

**4 - INDICATIONS AND CONTRAINDICATIONS**

**1) Indications**

- Post-op. pain therapy for a **maximum of 48 hours**.
- Biological safety risk due to analgesia if the duration of use is exceeded.
- Possible local and/or toxic effects on patient.

**2) Contre-indications**

- Lack of consent by the patient,
- Associated risks to the pain therapy.
- Proximity infection.
- Children less than 10 Kg.
- Use in epidural anesthesia.
- Known neurological abnormalities.

**5 - ADVERSE SIDE EFFECTS AND RESIDUAL RISKS**

**1) Adverse side effects**

- Lesion or laceration of tissues or vessels;
- Hematoma reactions.
- Air embolism aponevrotic.

**2) Residual risks**

- Risk of reduction or no flow due to incorrect placement or connection of the catheter;
- Risk of infection by rupture or forgetting of the catheter in situ at withdrawal;
- Local and/or toxic risks to drugs by overdose.
- Risk of chondrolysis.

**6 - STERILIZATION**

The **Profilis PaincAth®** set is supplied sterile in double packaging.

Sterilization by Ethylene Oxide.

Do not re-sterilize

Resterilization may result, in a failure of package integrity, a risk of cross-infection, a failure of performance of the device, a degradation of the material and cause endotoxic reactions due to residues.

The device contamination can lead to lesions, disease, up to the patient death. Do not re-sterilize the catheter whatever the process.

**The manufacturer cannot be liable of any damage coming from the re-use of the medical device or any of its accessories.**

The content of the conditioning is sterile as long as it is neither opened nor damaged.

Do not use if conditioning is damaged.

Infectious risk.

Expiry date: 4 years and 11 months.

**6 - INSTRUCTIONS FOR IMPLANTATION**

**Only experimented people can set up this device.**

Never set up this device by the intravenous way for the diffusion of analgesia. **Death risk.**


In orthopedic surgery and traumatology, respect the position of the catheter in the relation to the bone. **Risk of chondrolysis**

Two implantation techniques for the catheter **Profilis PaincAth®**:

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1. From the surgical cavity to the outside. Connect the catheter on the tunneler.

2. Pull the tunneler with its catheter outside the surgical cavity up to the marking of the catheter: position of the last perforation
3. Connect and secure the removable Luer connector on the catheter.

Do not screw excessively the Luer connections. Fissure risk.

4. Rinse the catheter with a sterile physiological solution using a syringe.

**b) Technique with peel-away introducer needle (R, QR versions) :**

- 1 From the outside of the surgical cavity.
- 2 Set up in the opened surgical cavity **ALL** the perforation length (closed distal end up to the marking of the catheter: position of the last perforation).
- 3 Rinse the catheter with a sterile physiological solution using a syringe.

For the 2 techniques of implantation:

For rinsing, never use a syringe inferior to 10 CC. Risk of burst of the catheter by overpressure.

**a) Close the surgical cavity.**

**b) Fix the catheter using the sterile fixing device provided.**

Do not use suture for fixing the catheter.

Risk of catheter pinch-off.

**c) Connect the catheter to the perfusion line.**

Always take care of the asepsis of the injection site. Do not screw excessively the Luer connections. Fissure risk.

**7 - ELIMINATION OF DEVICE AND ACCESSORIES**

At the end of treatment, the catheter will be softly removed by the authorized staff.

The physician will make sure that the catheter has been entirely removed (distal end marked). The discard of the device and its accessories must follow the protocol of the discard of medical wastes DASRI (Infections Risks Care Activities Waste) related to the centers, and must not endanger the safety or health of the patient the users or any other person until complete destruction. Infectious risk by accident of exposure to blood caused by puncture or cutting accessory.

In case of inefficiency of pain treatment with our Profilis PaincAth catheter and the analgesic, a substitution treatment chosen by the medical profession will have to be performed on the patient. The catheter will be removed. Any serious incident involving the device will be notified to the manufacturer and the competent authority of the member state in which the user and / or patient is established;

Kit unitaire / Unitary set		Kit par 10 / Box of 10 sets	
P500-02.5 S/R/P		P500-02.5 QS/QR/QP	
P500-05 S/R/P	P1000-05 S/R	P500-05 QS/QR/QP	P1000-05 QS/QR
P500-07.5 S/R/P		P500-07.5 QS/QR/QP	
P500-10 S/R	P1000-10 S/R	P500-10 QS/QR	P1000-10 QS/QR
P500-15 S/R	P1000-15 S/R	P500-15 QS/QR	P1000-15 QS/QR
P500-20 S/R	P1000-20 S/R	P500-20 QS/QR	P1000-20 QS/QR
P500-25 S/R	P1000-25 S/R	P500-25 QS/QR	P1000-25 QS/QR
P500-30 S/R	P1000-30 S/R	P500-30 QS/QR	P1000-30 QS/QR

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- Place des infiltrations en péroratoires. Injection unique vs cathétères vs systèmes à libérations prolongées. Marc Beaussier, Mourad Aïssou. Département d'Anesthésie-Réanimation chirurgicale, Hôpital St-Antoine. Groupe Hospitalier Est-Parisien. APHP, Paris - Sorbonne Universités. UPMC Univ Paris 06. France.Le Congrès. Medecins. Conférence d'Actualisation. © 2014 SFAR.
sont des marques déposées par / *are trade marks by* :
est une marque déposée de / is a trademark of / es una marca depositada por /
est une marque déposée de / is a trademark of