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
INTRODUCTION

The implantable catheter ports DistricAthe® are vascular access devices.

The wide range of the DistricAthe® implantable catheter port allows the doctor a choice adapted to the patient's morphology related to his future treatment.

The access way is percutaneous or surgical. The most used location is thoracic, the catheter path being a vein (by example : cephalic, sub-clavian, jugular...).

- Please read our technical brochure Doc.DistricAthe® to choose the best device for the patient in line with the patient's treatment and morphology.

This implanted medical device consists of an implantable port, a catheter and a connection ring.

DistrictAthe® ports are available in solid Titanium or Titanium / Polysulfone. Whatever the model, all the surfaces in contact with the injected products are made of Titanium. These injections are made by a non-coring needle (for example Huber tip).

TARGET POPULATION : Patient weighing over 10 kg and age over 2 years.

DEVICE PERFORMANCES AND BENEFICES

► Clinical performances :

- The DistricAthe® implantable catheter chamber allows repeated and prolonged direct access to the intravenous vascular system;

► Clinical Benefits :

- It makes it possible to conserve the venous capital of the patient ;
- The puncture in the port through the skin is less painful than a puncture in the vein ;
- Completely implanted under the skin, the first septic bulwark, the DistricAthe® implantable catheter port allows the patient a normal social life (shower, bathing, non-bath, physical exertion, ...) by minimizing the risk of infection compared to an exteriorized implanted system.

► Mechanical Performances :

- The titanium material combines resistance, lightness, biocompatibility and inertia with therapeutic products (Standard for surgical implants ISO 5832-2 Unalloyed titanium).
- The silicone septum ensures a tightness of 200 kPa (Standard ISO 10555-6 § 4.5.3 Absence of leakage, for 1000 punctures per cm²) with a Huber needle of 22G, that is : 1300 punctures for the standard range / 800 punctures for the small model range / 600 punctures for the junior range.
- The catheters are available in Silicone or Polyurethane. They can be pre-connected. The start of the marking is 10cm on the proximal side. Marking encrypted every 5cm, separated by dots every cm up to 45cm.
- The simple connection is made using a titanium fixing ring. The connection between chamber and catheter ensures a minimum hold of 5 Newton.
- The DistricAthe® port meets the most stringent pressure resistance requirements. The burst pressure of the system is oversized compared to the peak injection pressure without pressure.
- All the implanted elements are ORX and NMR compatible (up to 3 Tesla).

The summary of safety and performance characteristics is available on the EUDAMED website: <https://ec.europa.eu/tools/eudamed>

INDICATIONS AND CONTRA-INDICATIONS

1) INDICATIONS

The DistricAthe® implantable port in venous way is used for the administration of chemotherapy, antibiotic therapy and parenteral nutrition, as well as for the collection or injection of blood products.

- No other substance can be delivered to the DistricAthe® implantable port without medical validation. Uncontrolled biological risk for the patient.

2) CONTRAINDICATIONS

The exclusion criteria for the implantation of the catheter implantable port are :

Absolute contraindications :

- Patient's non-consent
- Infected areas, presence of bacteraemia, septicemia and peritonitis.
- Allergy to the materials making up the port and the catheter.
- The risks associated with local or general anesthesia.

Momentary contraindications:

- The anterior venous first placement by subclavian route.
- Upper cave syndrome.

Relative contraindications:

- Previously irradiated areas.
- Skin metastases.
- Mediastinal tumors.
- Disruption of blood clotting.

ADVERSE SIDE EFFECTS & RESIDUAL RISKS

1) ADVERSE SIDE EFFECTS

Percutaneous technique :

- Hematoma at the puncture site

Surgical technique :

- Hematoma in the subcutaneous pocket

2) RESIDUAL RISKS

- Infectious risk.
- Risk of extravasation of the drug.
- Obstruction of the catheter.
- Pneumothorax
- Rupture and migration of the catheter.
- Vein thrombosis.
- Skin ulceration and necrosis with regard to the implantable port.
- Respiratory distress (hydrothorax).

READ CAREFULLY BEFORE IMPLANTATION

SET COMPOSITION

The DistricAth[®] set is delivered sterile for single use with its accessories.

Do not re-use

The reuse exposes patient to infection risk.

Implanted parts :

- 1 Implantable port DistricAth[®] unique laser-marked SN number
- 1 Silicone Catheter or in polyurethane suited to the port exit
- 1 Connection ring suited to the catheter diameter

Accessories that can be integrated into the set (1) :

- 1 Puncture needle
18 Gauges - L = 70mm
or
1 Short puncture needle
21 Gauges - L = 40 mm

- 1 Huber straight needle
19, 20 or 22 Gauges - L = 20 to 35 mm
- 1 90° curved Huber needle
19, 20 or 22 Gauges - L = 20 to 35

- 1 Metallic guidewire
0.018" to 0.035"
L = 50 to 80 cm

- 1 Tuner
Ø 1.5 to 3 mm - L = 210mm

- 1 Peel-away Introducer standard or short
IPDM standard :
L = 205 mm ; Size = 5 to 10 French
IPDM short :
L = 100 mm ; Size = 5 to 7 French
IPG standard :
L = 215 mm ; Size = 7 and 9 French

- 1 Vein pick Length
L = 50 mm
Tip at 90° of 13mm

- 1 Scalpel
L = 14 mm - L = 138 mm

- 1 Syringe 10 ml centered Luer
L = 113 mm

- 1 Syringe 10 ml Luer Lock
L = 113 mm

- 1 Catheter Rinse - L = 34,5 mm

Sterile accessories under their own EC markings with their instructions for use than can be supplied on request :

- 190° Huber extension set
simple or with Y
19, 20 or 22 gauges - L=15/20/25mm

- 1 Introducer needle
12 Gauges - L = 80 mm

- 1 Cable for ECG connection ECG
L = 1 meter

STERILIZATION

The DistricAth[®] set is delivered sterile in double packaging

Ethylene oxide sterilized

Do not re-sterilize

Re-sterilization can compromise the integrity of the packaging, the risk of cross-infection, device performance failure, degradation of the material and lead to toxic reactions due to residues, or even death of the patient. Do not re-sterilize the port regardless of the process.

Expiration date Date : 4 years and 11 months
Do not use this device after the expiration date written on the package.

Do not use if the packaging is damaged. The contents of the package are sterile, as long as the package is neither usually open nor damaged, or does not show any weld breaks in the packages.

Contamination of the device can cause injury, illness, or even an infectious risk.

IMPLANTATION INSTRUCTIONS

1) PRECAUTIONS AND BAD REASONABLY PREDICTABLE LOCATIONS

Only qualified doctors can set up the DistricAth[®] implantable ports.

Before using this device, the doctor must have read this instruction for use in its entirety.

- The implantation of implantable catheter port requires surgical asepsis and suitable equipment :

Infectious risk for the patient

- Wrong orientation of the bevel of the puncture needle:

Risk of deformation of the guidewire spring (for example, in the case of a subclavian vein with a pronounced curvature). A blue mark on the base allows the orientation of the bevel.

- Blocking of the guidewire in the puncture needle :

When catheterizing a vein (Seldinger technique), never leave the guidewire immobilized in the puncture needle in order to avoid any risk of blockage of the latter by fibrin formation.

- Forced insertion of the introducer :

For some patients, to facilitate the penetration of the introducer, provide a slight skin incision near the guidewire to avoid the risk of kinking the introducer sheath.

- Weakening of the catheter by assembly using pliers/forceps :

Catheterization with risk of rupture over time and migration into the vascular system.

- Incorrect positioning of the connection ring :

Risk over time of catheter disconnection with migration of the catheter into the vascular network

- Rupture of the catheter by costo-clavicular pinching :

Never implant a subclavian catheter when there is a risk of costo-clavicular pinching. The risk is rupture of the catheter and its migration into the vascular network.

- Incorrect positioning of the catheter end:

Risk of device malfunction (no venous return by thrombosis of the catheter).

- Incorrect positioning of the port :

Shallow lodge. Risk of necrosis of the patient's skin
Lodge too deep. Risk of incorrect location by palpation.
Oversized lodge. Risk of overturning the port

2) IMPLANTATION SITE

Make sure you have all the necessary equipment for the implantation.

Before implantation, purge the port and the catheter.

The choice of the site where the implantable port will be placed may vary depending on the patient. You must select a location where the port rests on a resistant tissue (large muscle, bone surface).

The site must be adapted to the size of the port. Do not put the port too superficially to avoid any risk of necrosis of the skin. The site must be large enough to accommodate the entire port so that the scar is outside the puncture area.

The port has holes for the possible passage of stitches to secure it on the fascia.

CATHETER IMPLANTATION

The catheter can be placed either by surgical technique or by percutaneous technique. Depending on the technique, the catheter can take different access ways.

When using a percutaneous introducer, it must be adapted to the diameter of the catheter.

Two options are recommended for catheter placement:

- Leave the catheter connected to a syringe filled with physiological saline.
- Purge the catheter and clamp it at the distal end before cutting and inserting it into the vein, in order to avoid any blood rise (reduction of the risk of thrombosis).

In Europe, intravenous implantation accounts for 95% of implantable port cases.

To ensure the correct positioning of the proximal end of the catheter, the latter must be at the border of the superior vena cava, atrium, controlled by an image intensifier. As part of an implantation using ECG control technique, refer to our Noti033 supplementary notice.

Once the catheter is in place in the vein, provide an additional length of catheter to make the connection with the port.

The length of the catheter should be as short as possible in order to avoid the risk of kinking of the catheter at the exit of the port.

3) CONNECTION PORT / CATHETER

See description of the connection in page 2.

READ CAREFULLY BEFORE IMPLANTATION

4) OPERATING CONSEQUENCES

The follow-up operations are generally simple.
It must be ensured that the following local complications are not present:
Hematoma next to the room, it will subside in a few days: the room should not be used during this period if possible.
Infectious or inflammatory complications (redness, edema, pain...). The implantable catheter chamber should not be used until the skin is healthy again.
Skin sutures will be removed between the 7th and 10th day.

⚠ If water-absorbable threads are used, the scar should not be wetted for 7-10 days. Risk of wound opening.

The implantable system can be used immediately after placement if the skin sutures are not directly located opposite the port.
A front chest X-ray should be prescribed. It will allow you to check the position of the catheter, the end of which should go down to the junction of the superior vena cava and the right atrium.
Fill in the card holder and the patient monitoring book

USAGE

1) PRECAUTIONS AND REASONABLY PREDICTABLE MISUSE

Use and maintenance may only be carried out by persons duly authorized and having received specific training.
Before using this device, the practitioner must have read this user manual in its entirety.

• The use of implantable catheter chambers requires site asepsis on the patient and suitable equipment:

⚠ Infectious risk for the patient

• Use of an unsuitable needle (coring for the septum):

⚠ It is mandatory to use straight or curved Huber needles. The use of conventional needles would degrade the septum by coring. Risk of catheter thrombosis due to perforation.

• Incorrect needle puncture outside the septum:

⚠ Risk of perforation of the catheter with thrombosis of the catheter. Impossibility of venous return.

• Use of an unsuitable syringe less than 10cc:

⚠ It is mandatory to use syringes with a volume of 10cc or more. The pressure exerted during an injection is greater the smaller the area of the syringe plunger.
The use of syringes of volume less than 10cc could generate overpressure a risk of bursting or rupture of the catheter with migration into the vascular system.

⚠ The use of syringes of volume less than 10cc could generate overpressure a risk of bursting or rupture of the catheter with migration into the vascular system.

• Incorrect catheter rinsing:

⚠ Risk of blockage of the catheter requiring anticoagulant therapy to make the system permeable or removal of the system if treatment fails.

• MRI examination with metal accessory:

⚠ Do not use incompatible metal accessories for NMR imaging. Risk of movement of the accessory or overheating on the patient.

2) INJECTIONS AND SAMPLES

The use of a protocol takes into account the following different steps:

- Injection of medicinal products.
 - System maintenance.
 - Prevention of complications related to the use of these systems (infection, rupture of the catheter, thrombosis).
- These recommendations are not set in stone. Depending on the publications, adaptations could be made in particular in the area of maintenance of the implanted device (heparinized serum or not and frequency of rinsing).

Asepsis the site before any surgical gesture, skin preparation of the patient is necessary. It is carried out with an antiseptic soap from the same family as the antiseptic used in the operating room.

Access to the port must be in very specific aseptic conditions:

- The operator must wear surgical clothing, a cap, a bib and sterile gloves.
- A team of practitioners trained in the use of implantable sites. Any injection or withdrawal is done following steps below:

- Location of the port: locate the septum of the room by palpation.
- Stitching of the Huber needle perpendicular to the skin. Crossing the septum offers resistance.

The needle should be pushed in until it stops at the bottom of the port. The correct positioning of the needle is checked by blood reflux immediately followed by rinsing by pulsed injection of a minimum of 10 ml of physiological saline.

⚠ The length of the needle must be adapted to the depth of the implantation of the chamber, and its diameter with the type of treatment.
• For blood samples, the first ml of blood collected must be discarded.

⚠ Withdraw the needle while continuing to inject with the syringe to maintain positive pressure. Reduced risk of catheter thrombosis.

⚠ When withdrawing the needle inserted into the septum, prevent a rebound effect. Risk of accidental stinging.

The use of a three-way stopcock reduces the risk of air embolism when changing syringes.

Before any injection, infusion or blood sample, the system must be rinsed with a physiological saline solution (greater than 10 ml).
Rinsing at the end of use makes it possible to make a lock of physiological serum or heparinized serum according to the medical protocol.

3) INJECTIONS

For continuous infusions, the septum puncture is usually performed with a curved Huber needle, connected to an extension tube and its 3-way stopcock. The fixation is ensured by an adhesive bandage. The occlusive dressing should allow access to the three-way stopcock.

In the case of continuous infusions over several days, the needles are generally changed on average every 48 hours. During these infusions it is important to respect the established protocol.

4) SYSTEM MAINTENANCE

The purpose of maintenance is to prevent the risk of system obstruction.

The flushing rate of the system is variable according to the protocols of the establishments. Generally, the system is flushed with a physiological saline solution at implantation and then with a frequency of six to eight weeks in the absence of treatment.

a - Rinsing volume (ISO Standard 10555-6 54.5.4) physiological saline lock heparinized serum lock (See table p.23)

It will be essential to carry out a pulsed injection of at least 10 ml of physiological saline, followed by a lock of 5 ml of physiological saline or heparinized serum in slow injection.

NB : Doctors J.L. THEARD and S. ROBARD from the René Gauducheau Center in NANTES *, showed that rinsing with non-heparinized physiological saline does not increase the frequency of complications of injection site or catheter thrombosis. However, it eliminates the risk of overdose, product error and heparin intolerance. Finally, it reduces the cost and care of home maintenance by nurses, as well as the risk of infection by manipulations.

* J.L. THEARD, S. ROBARD. «Evolution of the protocols for using implantable chambers» Ann Fr Anesth Réanim, 1995; 14: 535

Between these infusion sessions, the patient continues his activities without constraints, except those advised against by his attending physician. Possible discomforts can be those caused by the wearing of seat belts and improper exercise.

b - Special cases

• In any case, never inject medicinal products (cytotoxic) without having checked the permeability of the system.

- If an injection is impossible, proceed to a radio control to find out the cause :
- If the catheter is thrombosed: use of anticoagulants (technique to be used only by qualified doctors).
- If the catheter is cut: removal and replacement by qualified doctors.

IMPLANTATION DURATION AND WITHDRAWAL

The duration of implantation of the DistricA[®] implantable catheter port depends on the duration of treatment.

- It is however recommended to remove the device according to several criteria :
- After the end of treatment depending on the psychological state of the patient.
- In the event of foreseeable non-treatment of more than 6 months following the establishment's protocol.

In addition, removal of an implantable port may be necessary, although the implant is not involved, due to side effects, such as local infection, tissue necrosis or even venous thrombosis.

The removal of this medical device must be done according to the rules of surgical asepsis with the appropriate equipment.

⚠ When removing the implantable port, remember to remove all the elements. Make a secure withdrawal of the catheter to avoid rupture of the catheter. Risk of migration.

⚠ The disposal of the device and of the waste must be done according to the protocol for the elimination of medical waste DASRI (Waste from Care Activities at Infectious Risks) linked to care establishments, and must not endanger safety or health, of the patient, the users, or any other person until its complete destruction.

⚠ Any serious incident occurring in connection with the device will be the subject of a notification to the manufacturer and to the competent authority of the Member State in which the user or /or patient is established.

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est une marque déposée de / is a trademark of / es una marca depositada por / é uma marca registrada de