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
READ CAREFULLY BEFORE IMPLANTATION

INTRODUCTION

Implantable catheter ports DistricAth® are vascular access devices. These devices allow the administration of medications and the withdrawal or the injection of blood products. These injections are performed using a non-carving needle (by example Huber tip).

DESCRIPTION AND DEVICE PERFORMANCE (ISO 10555-6 STANDARD)

This medical device is composed of an implantable port, a catheter and a connection ring. Implantable catheter ports DistricAth® are available in solid Titanium or Titanium/Polysulfone (Ti/PSU). **Whatever the reference, all areas in contact with the injected products are in solid Titanium.** The material is combining resistance, lightness, biocompatibility and inertia to therapeutic products.

Catheters are available in Silicone or Polyurethane. Catheters can be preconnected. The number 0 of the marking is at the proximal end.


Implantable catheter ports DistricAth® are allowing a repeated and extended access to the vascular system (intravenous or intra-arterial).

The silicone septum ensures a tightness of 200 KPa (ISO 10555-6 54.5.3). No leak for 1000 punctures by cm² with a Huber needle of 22G.

The simple setting is performed by a securing titanium ring. The DistricAth® port does answer to the strictest requirements regarding pressure resistance. The burst pressure of the system is oversized compared to the peak pressure, with no pressure of injection. The port-catheter connection does ensure a minimal holding of 5 Newton. Implantable catheter ports range DistricAth® is NMR (Nuclear Magnetic Resonance) and X-RAY compatible (up to 3 Teslas). All catheters are radio opaque (X-RAY). They are graduated from 10cm to 45cm, every 5cm separated by a point every cm.

SET COMPOSITION

The DistricAth® set is provided sterile, disposable with all its accessories.

Implanted parts :

- 1 DistricAth® implantable port, n° SN laser engraved
- 1 Silicone or polyurethane catheter suited to the port exit
- 1 Connection ring suited to the port exit

Accessories that could be in the Set (1) :

- 1 Puncture needle
18 Gauges - 70 mm

- or short Puncture Needle
21 Gauges - 40 mm

- 1 Huber straight needle
19, 20 or 22 Gauges
Length = 20 to 35 mm

- 1 Huber curved needle
19, 20 or 22 Gauges
Length = 20 to 35 mm

- 1 Guidewire
Length = 50 to 80 cm
Ø 0.018" to 0.035"

- 1 Peel-away Introducer standard or short
IPDM standard :
L = 205 mm ; Size = 5 to 10 French
IPDM short :
L = 100 mm ; Size = 5 to 7 French
IPG standard :
L = 215 mm ; Size = 7 and 9 French

- 1 Tunneler
Length = 210 to 218 mm
Ø 1,5 to 3 mm

- 1 Vein Pick
Length = 50 mm
Tip at 90° of 13 mm

- 1 Scalpel
Length = 138 mm
Wight = 14 mm

- 1 Syringe 10 ml
Luer Centred
Length = 113 mm

- 1 Syringe 10 ml
Luer Lock
Length = 113 mm

- 1 Syringe 10 ml
Luer decentred
Length = 113 mm

- 1 ECG Cable
Length = 1 meter

- 1 Rinse Catheter
Length = 34,5 mm

- 1 2nd Connection Ring

- 1 2nd Huber straight or curved needle

- 1 2nd syringe

- 1 2nd Guidewire

Sterile accessories with their EC marking provided upon request :

- 1 Huber Extension Set 90° simple or with Y connection

- 1 Peel-away Introducer Needle Picture

(1) Technical specifications file available for each accessory upon request

INDICATIONS AND CONTRAINDICATIONS

1) INDICATIONS

The access way is percutaneous or surgical. The main port implantation area is thoracic, with the catheter introduced into a vein (i.e. cephalic vein, jugular vein, subclavian vein...).

Kindly read our Technical file Doc.districtAth® to choose the appropriate device for the chosen treatment and the patient morphology.

The implantable port DistricAth® by venous way, is used for the administration of medications, chemotherapy, antibiotherapy and parenteral nutrition as well as the withdrawal or the blood products injections.

Check for physico-chemical incompatibilities that could exist between several medications solutions injected.

2) CONTRAINDICATIONS

The exclusion criteria for placing an implantable catheter port are :

Absolute contra-indications :

- The patient disagreement.
- Infected areas, presence of bacteria in blood, septicaemia and peritonitis.
- Allergy to port and catheter materials.
- Related risks to local or general anesthesia.
- For children of less than 10KG and less than 2 years old.

Temporary contra-indications :

- Earlier implantation by vein access through subclavian way.
- Superior cava syndrome.

Relative contra-indications :

- Earlier irradiated areas.
- Cutaneous metastasis.
- Mediastinal tumours.
- Blood coagulation troubleshooting.

3) POSSIBLE COMPLICATIONS

Percutaneous technique :

- Hematoma at the Pneumothorax puncture.

Surgical technique :

- Hematoma in the subcutaneous pocket.

4) SIDE EFFECTS

- Risk of infection.
- Risk of extravasation of the drug.
- Catheter obstruction.
- Venous thrombosis.
- Ulceration and skin necrosis above the implantable port.
- Respiratory distress (hydrothorax).

STERILIZATION

The DistricAth® set is delivered sterile and in double packaging.

Sterilized by Ethylene Oxide

The re-use of disposable devices exposes the patient or the user to infection risks. The device contamination can lead to lesions, disease, up to the patient death. Do not re-sterilize the implantable port whatever the process.

Do not re-sterilize

Deadline for use: 4 years and 11 months

Do not use this material after the expiry date which is stated on the carton box

The re-use will compromise the device integrity or will lead to the device default. **The manufacturer cannot be liable of any damaged coming from the re-use of the medical device or any of its accessories.**

The content of the conditioning is sterile as long as it is neither opened nor damaged.

Do not use if conditioning is damaged

READ CAREFULLY BEFORE IMPLANTATION

⚠ Do not use this device after the expiry date printed on the conditioning.

IMPLANTATION INSTRUCTIONS

1) PRECAUTIONS BEFORE USE

Only qualified physicians are authorized to implant DistricAth® ports.

⚠ When a catheter is implanted into the subclavian vein, it is forbidden to position the catheter into the costo-clavicular forcepts.

Be sure that you have all equipment needed for im-plantation.

⚠ During the catheterization of a vein (Seldinger technique), never allow the guide wire to become immobilized in the puncture needle to avoid any risk of blockage of it by fibrin formation.

If you do use a percutaneous introducer, check its compatibility with the chosen catheter (check the box ID sticker).

⚠ For some patients, to facilitate the penetration of the introducer, plan for a slight incision in the skin near the guidewire.

Purge the port and its catheter before implantation.

⚠ It is recommended, for the implantation of the catheter, two options :
- to leave the catheter connected to a syringe filled with physiological serum,
- to clamp the catheter before cutting it and then inserting it into the vein, to avoid any blood flushing back into the catheter (decrease of the risk of thrombosis).

Fill the CARD HOLDER identification and the follow-up notebook by putting the provided adhesive tracking sticker for the device reference and the set number n°SN.

2) IMPLANTATION SITE

The port location may vary with the patient. Choose a location where the port can be affixed on a resistant basis (large muscle, bone area).

The precise location shall be deep enough to receive the port set, so that the scar does not overlap it. Be careful not to implant the port too deep into the chosen site in order to easily locate it by palpation; or too superficially not to damage the skin (+/- 5mm). The port does include some perforations for stitches to affix it to the aponeurosis.

3) CATHETER IMPLANTATION

The catheter implantation may be performed by surgery or by percutaneous methods. Catheters may be placed in several locations depending upon the intended treatment area.

In Europe, the intravenous way does represent 95% of port implants. The catheter end is placed inside the superior cava vena, at the right atria orifice level. Once the catheter is inside the vessel, plan to have an extra catheter length to make the connection to the port exit. The catheter length shall be as short as possible to obtain the best flow rate.

Intra-peritoneal implantation remains a specific method for peritoneal carcinomas. The port can be implanted inside the ante superior iliac spine if the patients present a superior cava vena thrombosis. In that case the catheter should be introduced in the iliac vein with its extremity reaching the inferior cava vein.

4) PORT/CATHETER CONNECTION

See Connection steps in page 2

⚠ Do not use any pincer to set up the catheter on the port exit (catheter damage risk).

5) POST-SURGERY EFFECTS

Post-surgery effects are generally simple.

Be aware to check for any potential local post-surgery difficulties, such as:

Hematoma around the port, which will disappear after a few days; if possible, do not use the port during that time.

Infectious or inflammatory difficulties (blotch, edema, pain...). Do not use the port until the skin becomes healthy again.

Cutaneous stitches will be removed between the 7th and the 10th day.

⚠ If water resorbing stitches are used, do not wet the scar during 7 to 10 days.

The port set may immediately be used after the implantation if the stitches are not directly located above the port.

A front pulmonary radiography must be prescribed. It will help to check the catheter positioning which extremity shall go down to the superior cava vena and to the right atria junction.

USAGE

1) PRECAUTIONS OF USE

Using implantable ports with catheter requires a general surgical asepsis and appropriate equipments.

⚠ You must use straight or curved Huber type needles.

The use of normal needles will strongly damage the septum with perforations.

⚠ You must use syringes with a volume of 10cc minimum or more.

The use of smaller capacity syringe would damage the implantable port, since they are causing, with a blocked catheter, a rupture and a migration of that catheter into the blood circulation (the pressure due to an injection is increasing as the syringe piston area is decreasing).

The use of a three ways stopcock reduces the gaseous embolism risk when changing syringes.

Before and after each injection, perfusion or blood sampling, the system must be rinsed with a physiological serum solution (about 10 ml).

When finished, a physiological or heparinated serum bolt is performed according to the medical protocol.

⚠ Retrieve the needle while injecting with the syringe to maintain a positive pressure.

⚠ Maintenance has to be performed by qualified medical professionals after a specific training.

⚠ Only experimented personnel shall use this medical device. Before any use of this device, the personnel shall have read this Information For User booklet in full.

The use of a protocol shall take into account the following several steps :

- Drugs products injection
- Device maintenance
- Prevention of complications related to the use of these systems (infection, catheter break, thrombosis)

The recommendations can change. Depending upon the researches, the adaptations that could be made, especially in the maintenance area of the device (use of heparin and rinsing frequency).

2) INJECTIONS AND SAMPLINGS

A cutaneous preparation of the patient is necessary before any surgical action. This preparation is done with the same antiseptic soap family than those used in operating theater.

The installation of vascular devices must be done according to specific aseptic conditions:

- The operator must wear a surgeon's gown, cap, bib and sterile gloves.
- A well-trained physicians team specialized into ports implantation.

Each injection or sampling should be done in three steps:

- Locating the port: find the port septum by palpation, then hold the port base with your left hand if you are right - handed and vice-versa.
- Introduce the Huber needle perpendicularly to the skin. The crossing of the septum gives the sensation of piercing gum through.

The needle must be pushed until it reaches the port bottom.

A good needle positioning is checked by the occurrence of blood flowing back, which must immediately be followed by an injection of about 10 ml of physiological serum.

- Keep the piston pressure when you remove the needle to avoid blood flowing back into the catheter.

NOTE: Some physicians recommend to hold the syringe with one hand and the needle with the other hand to avoid the rebound effect. Other physicians prefer to hold the syringe with one hand and hold strongly the skin around the implantable port with the other hand.

⚠

- The Huber needle length shall be suited to the port implantation depth and the port diameter to the treatment type.
- For blood samplings, the first collected ml must be thrown away.

3) PERFUSIONS

For continuous perfusions, the septum puncture is usually performed with a curved Huber needle connected to an extension set. To maintain the set in place, one can use folded compresses and adhesive ribbon. The occlusal dressing must allow the access to the three ways stopcock. In case of long-time perfusions over several days, the needles are usually changed every 48 hours.

4) SYSTEM MAINTENANCE

Maintenance is to avoid the occlusion risk of the system. The rinsing frequency may change depending upon the hospital protocols.

Usually, the rinsing of the system with physiological serum is done during implantation, and then, with a frequency of six to eight weeks when no treatment is performed.

a - Rinsing volume (ISO Norm 10555-6 §4.5.4) physiological serum bolt heparinated serum bolt (see sheet in p.23).

It is essential to do a rapid injection of 10 ml of physiological serum, followed by a bolt of 5 ml of physiological serum or heparinated serum done by a slow injection. N.B.: About the use of heparin solution, a study carried out by doctors J.L. THEARD and S. ROBARD of the Center René Gauducheau in NANTES*, showed that the rinsing with a physiological salt solution not heparinated does not increase the frequency of the complications of thrombosis of the injection site or the catheter. On the other hand, it eliminates the risk of overdose, error of products and intolerance to heparin.

Finally, it allows to decrease the cost and the nursing care at home, as well as the risks of infection by the manipulations.

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⚠

- In any case, do not leave the needle in place if there is no perfusion or injection.
- The patient must behaves normally with the possibility to take showers or baths. The only inconvenience could occur when using car security belts.

- b - Particular cases
- In any case, never inject drug products (cytotoxics) without having checked the permeability of the system.
- In case of catheter obstruction:
 - Try to inject heparinated serum.
 - If this doesn't work, use an anti-coagulating (Heparine) in the port to break up the catheter obstruction (this technique must only be used by well-trained doctors).

⚠ NEVER USE A SYRINGE INFERIOR TO 10CC

- The replacement of the system can be considered if all these actions fail to work. In any doubt regarding the good positioning, a radiography will confirm the proper functioning of the system.

READ CAREFULLY BEFORE IMPLANTATION

IMPLANTATION LENGTH AND RETRIEVAL

The catheter implantable port DistricAth® implantation length is subjected to the treatment length and its efficiency.

However, it is recommended to retrieve the device depending upon several criteria :

- After a successful treatment: a treatment can be up to 3 months (chemotherapy) and until 5 years (viral treatment).
- If no treatment is planned for at least 6 months.

Moreover, the implantable port retrieval can be necessary, even if the device is not the issue, due to secondary effects, such as a local infection, a tissue necrosis or a venous thrombosis.

During the retrieval process, check to retrieve all implanted components. Proceed with a secure catheter retrieval to avoid any break of this one. Migration risk.

The retrieval of this medical device shall be performed following all surgical asepsis rules with the appropriate equipment.

The removal of the device and waste must be done according to the ACIRW medical waste removal protocol (Activities of Cares with Infectious Risks Waste) related to the centers, and must not endanger the safety or health of patients, users, or anyone else until its complete destruction.

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est une marque déposée de / is a trademark of / es una marca depositada por / é uma marca registrada de