districlass médical SA

endogAst[®]

PROTHÈSE INTRA GASTRIQUE AVEC
INTRA GASTRIC PROSTHESIS WITH
PROTESI INTRAGASTRICA CON
PRÓTESIS INTRAGÁSTRICA CON

districAth[®]

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TABLE OF CONTENTS

1- INTRODUCTION	Page 09
2- CAUTIONS	Page 09
3- SET PRESENTATION	Page 09
4- WARNINGS	Page 11
5- ELIGIBILITY CRITERIA	Page 12
6- ATIIP TECHNIQUE Short description	Page 12
7- CLINICAL PLAN PROPOSAL Summary	Page 14
8- FOLLOW-UP PLAN	Page 15
9- RETRIEVAL PROCESS	Page 15

1- INTRODUCTION

The intra gastric Prosthesis is an implantable medical device during 1 year for the treatment of patients with morbid or severe obesity for a temporary or permanent treatment.

The Prosthesis is implanted using a surgical and endoscopic procedure.

The prosthesis is progressively insufflated with air using a implantable Port, to reach it's nominal volume. The Prosthesis is interfering with the satiety control process and is inducing the reduction of meal introduction and weight loss.

2- CAUTIONS

Only qualified physicians trained to the ATIP Technique and duly licensed to practice in their country in accordance with the laws in force in this country can use and implant this Prosthesis.

Districlass médical SA will neither be responsible for any implantation consequences nor for a person who did not respect the above statement.

- Single use **STERILE** device.
- Ethylene oxide sterilization.
- **Check** the packaging before use.
- **Check** the expiry date on the sticker. Never use a device with a passed expiry date.
- **Do not clean. Under no circumstances, these products must be re-sterilized and reused** causing a major biological hazard for the patient such as bacteria, virus, infection ...
- Store in a cool, dark, dry place.
- Do respect all usual protocols in place in operating rooms to maintain the device sterility.
- Read carefully the user manual before implantation.

3- SET PRESENTATION

Packaging:

Medical device with its implantation accessories single use, **STERILE**.

Sterilization by Ethylene Oxide **EO**

Expiration: 5 years.

- 1 intra gastric prosthesis
 - Primary conditioning, in a single PETG blister thermosealed with medical grade cover (Tyvek®) and inserted into a PE pouch, thermosealed with medical grade cover (Tyvek®) or labelled unitary, in hard thermosealed double blister PETG with medical grade cover (Tyvek®).
 - Secondary conditioning, in labelled unitary box with 1 User manual, 1 card holder and tracking control stickers for the serial number SN.
- 1 sterile implantation tube conditioned under thermosealed peel away envelop.
- 1 sterile 60 ml syringe conditioned under thermosealed peel away envelop.

Features:

Endo 25D and Endo 40D versions

Weight:


- Prosthesis alone = 11 grams (Endo 25D)
12 grams (Endo 40D)
- Port with connection ring
- Compatible MRI device, radiopaque (X-RO)

Overall dimensions:

Oval Prosthesis

inflated with 250 ml (Endo 25D)

400 ml (Endo 40D)

The Set: (see picture below)

A - Intra gastric **endogAsti**[®] Prosthesis

- 01 1 Intra gastric Prosthesis
- 02 1 lasso wire in stainless steel
- 03 1 puncture cannula
- 04 1 3 ways Huber needle extension Set
- 05 1 tunneller
- 06 1 Implantable Port
- 07 1 simple stomach attach
- 07bis 1 double stomach attach
- 08 1 one way valve
- 09 1 connection ring
- 10 1 scalpel


B

- 11 1 implantation tube sterile with hook

C

- 12 1 sterile 60ml syringe

1 User Manual
4 tracking stickers
1 card holder

4- WARNINGS

The Prosthesis Implantation associated risks are similar to usual risks associated to any surgery or surgical acts, or anesthesia. Each patient may differently react to the surgery.

The patient shall have before an interview with the referent doctor, a nutritionist and eventually a psychiatrist.

Before the implantation, the patient shall sign a Patient Informed Consent with the implantation Doctor.

4.1 The setting error risk by the endoscope percutaneous technique is depending upon each patient. (Using this technique, a mini coelioscopy helps to prepare the anterior stomach face and allows the puncture needle to pass through with a maximum security level. Associated risks are the ones of a mini coelioscopy). It is under the sole responsibility of the implantation person to check that no other organ can be damaged at the time of the puncture; especially the liver hypertrophy common with obese people which has to be controlled with this type of patients. A scan is necessary to evaluate the liver dimensions.

4.2 Haemorrhagic risks at the puncture point.

There are no or little associated haemorrhagic risks with a puncture located at the left epigastria region. Each patient does differently react to the surgery: a minor hemorrhage may happen. If this phenomenon does happen, one must postpone the implantation of the **endogAst**[®].

4.3 Infectious risks.

An oozing gastric fluid outside the body or along the subcutaneous catheter may happen in particular cases and lead to a local skin inflammation. This kind of inflammation/infection is treated by local antibiotherapy. The possible encountered problems are identical to the ones of a percutaneous endoscopic gastrostomy.

To limit the infectious risks at the puncture point, an antibioprohylaxis is recommended the implantation day. It is imperative to follow the technical details as described in the ATIP Technique.

4.4 Nauseas and throw-ups.

Correctly positioned (corpus fundus), the intra gastric Prosthesis shall neither provoke nausea nor throw-up. Nevertheless the volume variation of the Prosthesis does

allow to solve this potential phenomenon.

Patients shall be warned against the tolerance variability from one Person to another.

4.5 Other risks.

endogAst[®] is an innovative device, the migration shall be nil because of the Prosthesis fixation to the stomach wall.

Districlass médical SA took all the necessary care from the design to the manufacture of this intra gastric Prosthesis. Nevertheless the use of implantation surgical accessories (needle, scalpel, pince ...) may damage the Prosthesis.

When in use, the Prosthesis may require an addition of air to compensate the temperature variation phenomenon between the insufflated air and the human body (37°C).

It is imperative to use a needle with Huber's bevel to proceed to the insufflation or the deflation of the prosthesis. Other type of needle would damage the septum of the implantable port and the prosthesis would empty of its air.

4.6 Failure and damages.

The Patient shall be duly informed that this technique does provide an important help to loose weight, knowing that a Prosthesis efficiency tolerance always is dependant upon the nutritional habits change after implantation. Failure is therefore possible.

When implanting the Prosthesis, the surgeon shall check via the endoscope, the Prosthesis integrity and its good positioning. **Districlass Médical SA** can be kept under no circumstances responsible for this failure.

In case of deflagration of the Prosthesis in time, **Districlass Médical SA** cannot be liable, because the Company does not have any control over the Patient's social life, sports life or nutrition.

No financial damage or any other type of damages can be asked to **Districlass Médical SA** for:

- Implantation failure, treatment failure, Prosthesis retrieval failure or prosthesis replacement failure.
- Error of manipulation or maintenance.
- Non respect of inclusion criteria and of warnings.
- Prosthesis intolerance.

Districlass médical SA responsibility is limited to the Prosthesis replacement in case of problem. The problem validation will be performed by the Device analysis within the Quality Department of the Company in the strict observance of the 93/42/EEC Directive requirements.

5- ELIGIBILITY CRITERIA

FOLLOWING THE PROPOSED TECHNIQUE

Inclusion criteria:

(To be performed within 30 days prior to implant).
Criteria are identical to the ones of obesity surgery.

Guidelines Consensus Conference - NIH - USA - 1991
for Obesity Surgery.

National Institute of Health consensus Development
conference statement 1991, March 25-27.
Am J Clin Nutr 1992 ; 55:615S-9S 1991 - updated in
2004 at the ASBS consensus conference.

- **BMI ≥ 40** if patients strongly desire substantial weight loss, because obesity severely impairs the quality of their life.
- **BMI 35-39** with a documented co-morbidity. Included in this category are patients with high risk co-morbid conditions such as life threatening cardiopulmonary problems (e.g. severe sleep apnea, Pickwickian syndrome, obesity related cardiomyopathy, or severe diabetes mellitus).
- Other possible indications for patients with BMI's between 35 and 40 include obesity-induced physical problems that are interfering with lifestyle (e.g. musculoskeletal or neurologic or body size problems precluding or severely interfering with employment, family function and ambulation).
- **Age ≥ 18.**
- Obesity **since more than 5 years.**
- At least 3 previous medical treatments.
- No endocrine pathology.
- No psychiatric pathology.
- No malignant tumor or severe pathology.
- Patient motivated and available to participate to the treatment and long term follow-up.

Additional indications about the inclusion criteria: Morbid obese patients with more than 60 years and super-obese patients (BMI > 50) are able to be enrolled in the survey.

ATIIP can be performed either with the intent of a definitive therapy, or with the intent to reduce a too high perioperative risk before other surgical procedures for the treatment of morbid obesity (over 1 year).

Exclusion criteria:

- Large iatal hernia (> 5 cm).
- Gastric ulcer.
- Previous oesophageal and / or gastric surgery.
- Severe chronic gastritis.
- Any contraindication to an upper digestive endoscopy.
- Previous abdominal surgery that makes the procedure complex.
- Severe infections or illnesses of the abdominal wall (skin, subcutaneous layer, abdominal fascia and muscle).
- Pregnant woman.

6- ATIIP TECHNIQUE

SHORT DESCRIPTION


The ATIIP procedure is performed using a surgical and endoscopic procedure in two subsequent phases:

- 1) The Prosthesis is implanted in the body-fundus area of the stomach. The prosthesis implantation is using the Percutaneous Gastrostomy Endoscopic - (PEG-pull technique - concepts).
- 2) The Prosthesis is connected to a subcutaneous totally implanted system: catheter Port. General anesthesia (recommended) or sedation is required.

Phase n. 1

- Position the Prosthesis **01** inside the implantation tube using the hook **11**. Sterile gel is useful to help the sliding process **11**.
- Position the Patient on a table suited to the person' weight and morphology in a anti Trendeleburg (30°) position.

- Gastroscopic control and gastric insufflations to bring the stomach into apposition to the abdominal wall.


- The abdominal entry site is located between the median lines, 0,5-1 cm below the xifoideous process and the left costal margin that correspond to the insertion of the 7th rib to the sternum. This entry site does correspond to the anterior wall of the corpus-fundus (Labbe triangle) area.

- A 1,5-2cm length transversal skin incision and a complete subcutaneous dissection until the abdominal wall fascia, is done. **Precisely locate the depression made with the doctor's finger with the endoscope. This does prove the absence of any organ at the puncture point.**


- Placement of a needle cannula 03 into the stomach.
- Leaving the cannula in place, pull out the needle. Very important ⚠
- Passage of the guidewire 02 into the stomach. Place the endoscope lasso around the cannula to facilitate the guidewire catching.

- Pull the guidewire out the patient's mouth while maintaining strongly the cannula in place.
- Make a buckle with the two wires, the one of the Prosthesis 01 and the one of the guidewire 02.

- Get the prosthesis down into the stomach by pulling the guidewire while pushing the tube 11 in the eosophagus. Always by maintaining the cannula in place.
- When passing the gastric wall, the catheter shall be pulled perpendicularly to the abdomen.


- Cut the catheter and throw away everything (except external 10-15 cm of catheter, to be able to perform the subcutaneous tunnellation).
- Endoscopic check of the prosthesis proper positioning inside the stomach.

Phase n. 2

- Change of the surgeon's gloves with a new pair ⚠
- A 3cm length slightly oblique skin incision is made about 10cm far from the median incision, 0.5 cm below the area of left rib cage. A subcutaneous pocket is performed for the port allocation 06 (remote from the scan) and the silicone cylinder 07 or 07 bis.
- The polyurethane catheter is tunnelled subcutaneously 05 from the median to the lateral incision.


- The stomach attach simple **07** or double **07 bis** (1,5 cm diameter) is inserted over the polyurethane catheter and positioned subcutaneously over the abdominal fascia.

- Insert the ring **09** on the catheter.
- Connect the catheter to the port. Insert the catheter on the port (about 5 mm).
- Holding the set together (ring and catheter) push the ring until it reaches the port **06**.

- Inflate (250ml / 400ml of air, depending volume) the Prosthesis by using the Huber needle via the Port to check its good positioning inside the stomach; then deflate it.

- Insert the port **06** in the pocket and fix over the fascia (4 non absorbable stitches).
- Subcutaneous pediatric drains are positioned in the surgical incisions and removed 48 hours later.

The definitive Prosthesis' air inflation is performed on the ward, 48 hours after the operation using a Huber needle.

- The Prosthesis' air inflation is started with 150ml (Endo 25D) or 300 ml (Endo 40D) of air, depending volume.

On the next visits, 50ml of air can be further added, according to the patient's tolerance, until nominal volume.

7- CLINICAL PLAN PROPOSAL

SUMMARY

Pre- operative evaluation:

- Clinical and nutritional status evaluation.
- Psychological evaluation.
- Liver and abdominal ultrasound.
- Oesophagogastroduodenoscopy.
- Oral microflora examination (swab).

Three days before operation:

- Chlorhexidine digluconate 0,2% : start with 15 milliliters (mL) as a mouth wash for 30 seconds two times a day.

The day before operation:

- Start low weight heparin once a day until 15 days after the discharge of the patients.

The operation day:

- Intravenous empirical antibiotic prophylaxis is started the day of the operation (Ceftazidime and Teicoplanine) until the 24th post- operative hours.

Post-operative course:

- **1st** post-operative day : - Nothing by mouth.
 - Hydratation i.v.
 - Omeprazol by i.v. 40mg/ day/single dose.
- **2nd** post-operative day : - Prosthesis insufflation.
 - Abdominal X-Ray.
 - Starts liquids by mouth (water).
 - Remove the drain.
 - Hospital discharge.
 - Start semi-liquid diet.

- **Starts solid diet from the 10th post-operative day, with a suited and controled diet.**

Suggestions:

- Careful control of the surgical incisions is mandatory for the first 15 days after the operation.
- Avoid any abdominal workout during the first 2 post-operative months.
- It is mandatory to have a post operative clinical follow-up that includes physical activity and a substantial modification of the nutrition behaviour of the patient.

8- FOLLOW-UP PLAN

The following evaluation intervals are recommended for the Patients' follow-up:

1st month (+/- 7 days):

Medical chart.

2nd month (+/- 7 days):

Medical chart.

3rd month (+/- 7 days):

Medical chart.

6th month (+/- 14 days):

Medical chart - Abdominal X-Ray

Gastroscopy.

9th month (+/- 14 days):

Medical chart.

12th month (+/- 14 days):


Medical chart - Abdominal X-Ray

Gastroscopy.

- **Abdominal X-Ray and gastroscopy can be performed any time if the patient does complain with abdominal pains or with no more efficacy of the therapy.**
- **All the clinical evaluations and blood test that the clinical investigator considers necessary for an appropriate follow-up should be performed, even though not requested by the clinical plan of the study.**

9- RETRIEVAL PROCESS

The Prosthesis retrieval is conducted under the same conditions as the implantation via the endoscope, but under local anaesthesia.

- Fully deflate the Prosthesis (syringe 12 and valve 08) and control by endoscope the total prosthesis deflation.
- Incise at the puncture point and at the Port location site 06.
- Retrieve the Port and the silicone cylinder out of the lateral incision.
- Cut the catheter.
- Using the endoscope grip, retrieve the Prosthesis (prolin for security).

In case of incomplete deflation of the prosthesis, it is recommended to proceed with a total deflation by puncture or cut using the endoscope.

CAUTION :

The removal of the device and waste must be done according to the ACIRW medical waste removal protocol (Activities of Cares with Infectious Risks Waste) related to the centers, and must not endanger the safety or health of patients, users, or anyone else until its complete destruction.

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est une marque déposée de / is a trademark of / è un marchio depositato di