districlass médical SA

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INTRODUCTION

A pneumothorax is the presence of air in the pleural cavity, outside the pulmonary parenchyma.

The main cause of the spontaneous pneumothorax is the rupture of the small sub pleural parenchyma bubbles.

The pleural abrasion brush **netcAth®** is a medical device of class **Ila** for the treatment of recurrent spontaneous pneumothorax by mechanical abrasion.

DESCRIPTION

The brush **netcAth®** is moulded as a single piece mixing polyethylene / polyolefine white semi-flexible, biocompatible and radio opaque (X-Ray), **Free latex**.

Two distinct parts:

- A 13mm length extremity and 7mm wide, slightly bended, allows an efficient handling, by an endoscopic flat nozzle.
- A cylindrical half sphere head of diameter 10mm is covered by a multiple of small flexible pins with a scraper effect, is forming the brush part on a 17mm length. The height of the pins is 2mm.

PERFORMANCES

- The abrasion pleural brush **netcAth®** is a consumable medical device which allows a mechanical irritation of the pleura during opened or endoscopic surgery, in order to obtain a pleural pulmonary seam line in case of pneumothorax.

- The height of the pins has been validated in order to collect the parietal pleural fragments during abrasion.

- The extremity of the brush with lateral edges avoids a lateral swing during the handling and allows an efficient handling with an endoscopic flat nozzle 5mm wide.

- Its maximum diameter of 10mm allows an easy introduction with the gripper in a 10.5mm trocar.

- It can resist to shear strains up to 100 Newton during the handling.

- The material is biocompatible radio opaque (X-Ray).

SET COMPOSITION

The set includes 4 pleural abrasion brushes **netcAth®**. The **netcAth®** set is provided sterile and single use.

SET COMPOSITION

The set includes 4 pleural abrasion brushes **netcAth®**. The **netcAth®** set is provided sterile and single use.

Under no circumstances these products have to be re-sterilized, reconditioned and re-used entailing a major biological risk for the patient such as bacteria, virus, infection...

INDICATIONS AND CONTRAINDICATIONS

5) INDICATIONS

The surgical treatment is generally the pneumothorax recurrence, and includes two aspects:

- The surgical resection of the causal source, generally bullous dystrophies.
- A parietal pleura irritation to obtain a seam line of the pulmonary parenchyma, by mechanical pleural abrasion.

There is no target population apart contraindications and factors that might compromise the success of the operation.

The scientific studies for the treatment of the spontaneous pneumothorax show patients from 14 to 68 years old.

6) CONTRAINDICATIONS

- Tobacco.
- Scuba diving
- Flight travel during the 6 weeks after surgical operation.
- Sporting activity during the month following the surgical operation.
- Children less than 14 years old.

7) POSSIBLE COMPLICATIONS

The encountered complications are the same as those of the mechanical pleural abrasions used for the pneumothorax.

8) SIDE EFFECTS

- Pain.
- Pneumothorax recurrence.
- Bronchial obstruction and pneumonia.
- Infection risk.
- Respiratory distress (hydrothorax or hemothorax)
- Blood clots forming into the veins (thrombosis).

STERILIZATION

The **netc[△]th[®]** is provided sterile in a double packaging. Sterile device.

STERILE EO Sterilization by Ethylene Oxide.

The use has to be made in the respect of the usual procedures of the operating room to protect the sterility.

The re-use of single use devices exposes the patient or the user to infection risks. The device contamination can lead to lesions, disease, up to the patient death. Do not re-sterilize the device whatever the process.

Do not re-sterilize.

The re-use will compromise the device integrity or will lead to the device default.

Do not use if conditioning is damaged.

Do not use after the expiry date mentioned on the conditioning.

Expiry date: 5 years.

USAGE

Before using the pleural abrasion brush **netc[△]th[®]**, the surgeon must be able to justify experience and training in the field of thoracic surgery. Only experimented surgeons shall use this device.

Before any use of this device, the surgeon has to read carefully and entirely this user manual.

The surgeon grasps the **netc[△]th[®]** brush and positions it on the endoscopic gripper. This brush is introduced into the thorax through a trocar of 10,5mm. Under endoscopic control, the abrasion is performed to obtain a pleural weeping. When the brush is saturated by the pleural fragments, the surgeon is replacing it and renews the process.


Easy to handle, thanks to its shape suited to the gripper, efficient during brushing due to its oval design and the multitude of pins, the pleural abrasion brush **netc[△]th[®]** is a simple solution, avoiding any recurrence during pneumothorax treatments.

Recommendation: a grooming by washing and/or aspiration is performed systematically after an irritation of the pleural cavity.

Advice (see drawings): according to the axis used, we do obtain a high flexibility at the **netc[△]th[®]** brush head (X axis), or a high hardness (Y axis). It does allow all abrasion possibilities, fine or strong.

X Axis

Y Axis

During the retrieval process, check that the brush has been entirely retrieved (Head, pins).

The remove of the device and waste must be done according to the ACIRW medical waste (Activities of Cares with Infectious Risks Waste) related to the centers, and must not endanger the safety or health of patients, or anyone else until its complete destruction.

RESPONSIBILITY

District medical SA, do certify that all the necessary care has been provided for the manufacturing of the device pleural abrasion brush **netc[△]th[®]**.

The manufacturer cannot be liable of any damages resulting from the re-use or re-sterilization of the medical device.

The content of the conditioning is sterile as long as it is neither opened nor damaged.

The manufacturer does not have any control over the selection of patients, the surgical handling, the conditions of use and the post-op. follow-up.

He cannot therefore guarantee no positive or negative effect.

The manufacturer disclaims all liability in case of damage, related directly or indirectly to the use of this device.

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est une marque déposée de / is a trademark of