districlass médical SA

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**1 - INTRODUCTION**

The multiperforated **Profilis PaincAth®** catheter is a medical device allowing the diffusion in situ of analgesic drugs for post-op. pain therapy. For the following molecules:

1. Bupivacaine (scientific references 1 ; 2 ; 3 ; 4)
2. Ropivacaine (scientific references 1 ; 2 ; 3 ; 4 ; 5 ; 6)
3. Lévocabupivacaine (scientific references 1 ; 2 ; 5)

**2 - DESCRIPTION AND PERFORMANCES OF THE DEVICE**

The multiperforated **Profilis PaincAth®** catheter is made of biocompatible radio opaque polyurethane.  
 Outside Ø = 1 mm ; Inside Ø = 0,5 mm ; Length = 50 cm or 100 cm.  
 The catheter is perforated with micro holes on 02.5 / 05 / 07.5 / 10 / 15 / 20 / 25 or 30cm, distal end side closed up to the 3 dashes  
 (See commercial references on a table at the end of the IFU).

**Diffusion area**

- The proximal extremity allows the connection to the perfusion line either by a pre-assembled connector or by a removable connector.
- Connection to an analgesic infusion line using its Luer Lock connector.
- Diffusion of the analgesia in situ **up to 48 hours maximum**

Biological safety risk:  
Possible local and/or toxic effects on patient.

- This kind of analgesia is an alternative to the use of morphine (side effects).
- Improvement of patient care post-op pain (Scientific references: 6; 7; 8).

- X-Ray and MR compatible 3 Teslas. **MR**

- Do not use the tunneler or introducer needle during NMR imaging.

- Free latex Free phtalate

**2 techniques of implantation:**

- \* References with S, QS or P, QP: tunneler technique
- \* References with R, QR: peel-away introducer needle technique

**3 - COMPOSITION OF THE SETS**

The **Profilis PaincAth®** set is supplied sterile, single use.

Do not re-use.

- The re-use of the single use device exposes the patient to infection risks.

**COMPOSITION OF THE S OR P SET**

- 1 X-Ray multiperforated catheter: OD = 1 mm ; ID = 0,5 mm ; Length : L = 50 cm or 100 cm.
- 1 Tunneler : S version (L = 210 mm) or P version (L = 100mm)

- 1 Removable Luer Lock connector

- 1 Luer Lock syringe of 10cc

**COMPOSITION OF THE R SET**

- 1 X-Ray multiperforated catheter: OD Ø = 1 mm ; ID Ø = 0,5 mm ; Length : L = 50 cm or 100 cm.

- 1 Peel-away introducer needle Ref. AIPM4 18 Gauges ; Length = 10 cm

- 1 sterile securement device, supplied under its own EC mark in each set. Ref. A004P 16/17G LOCKIT PLUS™

The references starting with Q (QS/QP/QR) mean that the box is containing 10 sets (see components of each set on the labels).

Possibility of a simultaneous diffusion on 2 catheters with the use of a Y Luer connector.

Accessory not supplied.

For safe administration, a 0,2µm antibacterial filter should be used on the infusion line.

Accessory not supplied.

**4 - INDICATIONS AND CONTRAINDICATIONS**

**1) Indications**

Post-op. pain therapy for a maximum of 48 h (See risk §2).

**2) Contre-indications**

- Lack of consent by the patient,
- Associated risks to the pain therapy.
- Proximity infection.
- Children less than 10 Kg.
- Use in epidural anesthesia.
- Known neurological abnormalities.

**3) Possible complications**

- Hematoma or allergic reactions.
- Air embolism aponevrotic.
- Lesion or laceration of tissues or vessels.
- Catheter pinch-off and reduced or absence of flow.
- Catheter break in situ during ablation.

**5 - STERILIZATION**

The **Profilis PaincAth®** set is supplied sterile in double packaging.

Sterilization by Ethylene Oxide.

Do not re-sterilize

Resterilization may result, in a failure of package integrity, a risk of cross-infection, a failure of performance of the device, a degradation of the material and cause endotoxic reactions due to residues.

The device contamination can lead to lesions, disease, up to the patient death. Do not re-sterilize the catheter whatever the process.

**The manufacturer cannot be liable of any damage coming from the re-use of the medical device or any of its accessories.**

The content of the conditioning is sterile as long as it is neither opened nor damaged.

Do not use if conditioning is damaged.

Infectious risk.

Expiry date: 4 years and 11 months.

**6 - INSTRUCTIONS FOR IMPLANTATION**

**Only experimented people can set up this device.**


Never set up this device by the intravenous way for the diffusion of analgesia. **Death risk.**

2 implantation techniques for the catheter **Profilis PaincAth®** :

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1. From the surgical cavity to the outside. Connect the catheter on the tunneler.

2. Pull the tunneler with its catheter outside the surgical cavity up to the marking of the catheter: position of the last perforation
3. Connect and secure the removable Luer connector on the catheter.

Do not screw excessively the Luer connections. Fissure risk.

4. Rinse the catheter with a sterile physiological solution using a syringe.

**b) Technique with peel-away introducer needle (R, QR versions) :**

- 1 From the outside of the surgical cavity.
- 2 Set up in the opened surgical cavity **ALL** the perforation length (closed distal end up to the marking of the catheter: position of the last perforation).
- 3 Rinse the catheter with a sterile physiological solution using a syringe.

For the 2 techniques of implantation:

For rinsing, never use a syringe inferior to 10 CC. Risk of burst of the catheter by overpressure.

- a) Close the surgical cavity.
- b) Fix the catheter using the sterile fixing device provided. Do not use suture for fixing the catheter.

Risk of catheter pinch-off.

- c) Connect the catheter to the perfusion line.

Always take care of the asepsis of the injection site. Do not screw excessively the Luer connections. Fissure risk.

**7 - ELIMINATION OF DEVICE AND ACCESSORIES**

At the end of treatment, the catheter will be softly removed by the authorized staff.

The physician will make sure that the catheter has been entirely removed (distal end marked). The discard of the device and its accessories must follow the protocol of the discard of medical wastes DASRI (Infections Risks Care Activities Waste) related to the centers, and must not endanger the safety or health of the patient the users or any other person until complete destruction. Infectious risk by accident of exposure to blood caused by puncture or cutting accessory.

In case of inefficiency of pain treatment with our Profilis PaincAth catheter and the analgesic, a substitution treatment chosen by the medical profession will have to be performed on the patient. The catheter will be removed

Any serious incident involving the device will be notified to the manufacturer and the competent authority of the member state in which the user and / or patient is established;

Kit unitaire / Unitary set		Kit par 10 / Box of 10 sets	
P500-02.5 S/R/P		P500-02.5 QS/QR/QP	
P500-05 S/R/P	P1000-05 S/R	P500-05 QS/QR/QP	P1000-05 QS/QR
P500-07.5 S/R/P		P500-07.5 QS/QR/QP	
P500-10 S/R	P1000-10 S/R	P500-10 QS/QR	P1000-10 QS/QR
P500-15 S/R	P1000-15 S/R	P500-15 QS/QR	P1000-15 QS/QR
P500-20 S/R	P1000-20 S/R	P500-20 QS/QR	P1000-20 QS/QR
P500-25 S/R	P1000-25 S/R	P500-25 QS/QR	P1000-25 QS/QR
P500-30 S/R	P1000-30 S/R	P500-30 QS/QR	P1000-30 QS/QR

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- Place des infiltrations en peropérateurs. Injection unique vs cathéters vs systèmes à libérations prolongées. Marc Beaussier, Mourad Aissou. Département d'Anesthésie-Réanimation chirurgicale, Hôpital St-Antoine. Groupe Hospitalier Est-Parisien. APHP. Paris - Sorbonne Universités. UPMC Univ Paris 06. France. Le Congrès. Médecins. Conférence d'Actualisation. © 2014 SFAR.
est une marque déposée de / is a trademark of