districlass médical SA

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1 - INTRODUCTION

The multiperforated **Profilis PaincAth®** catheter is a medical device allowing the diffusion in situ of analgesic drugs for the treatment of patient pain in post-op, dedicated to surgeries such as: rachis, pediatric, orthopedic, laparotomy, urology, visceral, gynecology, caesarian,

The multiperforated **Profilis PaincAth®** catheter is made of biocompatible radio opaque polyurethane, free latex and free phthalate

Outside Ø = 1 mm; Inside Ø = 0,5 mm; Length = 50 cm or 100 cm.
The catheter is perforated with micro holes on 2.5 / 5 / 7.5 / 10 / 15 / 20 / 25 or 30cm, distal end side closed up to the 3 dashes.

(See commercial references on a table at the end of the IFU).

TARGET POPULATION: Patient with a weight over 10 kg without age limit.

2 - PRODUCT PERFORMANCES AND BENEFICES

>Clinical Performances :

- The multiperforated catheter **Profilis PaincAth®** is a medical device allowing the in situ diffusion of analgesic drugs up to **48h00 maximum** for the treatment of post op pain for the following molecules:

1. Bupivacaine (scientific references 1 ; 2 ; 3 ; 4)
2. Ropivacaine (scientific references 1 ; 2 ; 3 ; 4 ; 5 ; 6)
3. Lévocabupivacaine (scientific references 1 ; 2 ; 5)

>Clinical Benefices :

- Continuous Analgesia effective for 48 hours (Scientific References: 6 ; 7 ; 8).
- Alternative to opioid consumption and their associated side effects
- Faster post-operative recovery;
- Decrease in the hospitalization length.

>Mechanical Performances :

- The proximal extremity allows the connection to the perfusion line either by a pre-assembled connector or by a removable connector.
- Connection to an analgesic infusion line using its Luer Lock connector.

• X-Ray and MRI compatible.

Do not use the tunneler or introducer needle during MR imaging.

3 - COMPOSITION OF THE SETS

The **Profilis PaincAth®** set is supplied sterile, single use.

Do not re-use.

The re-use of the single use device exposes the patient to infection risks.

2 techniques of implantation:

- * References with S, QS or P, QP: tunneler technique
- * References with R, QR: peel-away introducer needle technique

COMPOSITION OF THE S OR P SET

- 1 X-Ray multiperforated catheter: OD = 1 mm ; ID = 0,5 mm ; Length : L = 50 cm or 100 cm.
- 1 Tunneler : S version (L = 210 mm) or P version (L = 100mm)

- 1 Removable Luer Lock connector

- 1 Luer Lock syringe of 10cc

COMPOSITION OF THE R SET

- 1 X-Ray multiperforated catheter: OD Ø = 1 mm ; ID Ø = 0,5 mm ; Length : L = 50 cm or 100 cm.

- 1 Peel-away introducer needle
Ref. AIPM4
18 Gauges ; Length = 10 cm

- 1 sterile securement device, supplied under its own EC mark in each set. Ref. A004P 16/17G LOCKIT PLUS™

The references starting with Q (QS/QP/QR) mean that the box is containing 10 sets (see components of each set on the labels).

Possibility of a simultaneous diffusion on 2 catheters with the use of a Y Luer connector. Accessory not supplied.
For safe administration, a 0,2µm antibacterial filter should be used on the infusion line. Accessory not supplied.

4 - INDICATIONS AND CONTRAINDICATIONS

1) Indications

Post-op. pain therapy for a **maximum of 48 hours.**

Biological safety risk due to analgesia if the duration of use is exceeded.
 Possible local and/or toxic effects on patient.

2) Contre-indications

- Lack of consent by the patient,
- Associated risks to the pain therapy.
- Proximity infection.
- Children less than 10 Kg.
- Use in epidural anesthesia.
- Known neurological abnormalities.

5 - ADVERSE SIDE EFFECTS AND RESIDUAL RISKS

1) Adverse side effects

- Lesion or laceration of tissues or vessels;
- Hematoma reactions.
- Air embolism aponevrotic.

2) Residual risks

- Risk of reduction or no flow due to incorrect placement or connection of the catheter;
- Risk of infection by rupture or forgetting of the catheter in situ at withdrawal;
- Local and/or toxic risks to drugs by overdose.
- Risk of chondrolysis.

6 - STERILIZATION

The **Profilis PaincAth®** set is supplied sterile in double packaging.

Sterilization by Ethylene Oxide.

Do not re-sterilize

Resterilization may result, in a failure of package integrity, a risk of cross-infection, a failure of performance of the device, a degradation of the material and cause endotoxic reactions due to residues.

The device contamination can lead to lesions, disease, up to the patient death. Do not re-sterilize the catheter whatever the process.

The manufacturer cannot be liable of any damage coming from the re-use of the medical device or any of its accessories.

The content of the conditioning is sterile as long as it is neither opened nor damaged.

Do not use if conditioning is damaged.

Infectious risk.

Expiry date: 4 years and 11 months.

6 - INSTRUCTIONS FOR IMPLANTATION

Only experimented people can set up this device.


Never set up this device by the intravenous way for the diffusion of analgesia.

Death risk.
In orthopedic surgery and traumatology, respect the position of the catheter in the relation to the bone. **Risk of chondrolysis**

Two implantation techniques for the catheter **Profilis PaincAth®**:

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1. From the surgical cavity to the outside. Connect the catheter on the tunneler.

2. Pull the tunneler with its catheter outside the surgical cavity up to the marking of the catheter: position of the last perforation
3. Connect and secure the removable Luer connector on the catheter.

Do not screw excessively the Luer connections. Fissure risk.

4. Rinse the catheter with a sterile physiological solution using a syringe.

b) Technique with peel-away introducer needle (R, QR versions) :

- 1 From the outside of the surgical cavity.
- 2 Set up in the opened surgical cavity **ALL** the perforation length (closed distal end up to the marking of the catheter: position of the last perforation).
- 3 Rinse the catheter with a sterile physiological solution using a syringe.

For the 2 techniques of implantation:

For rinsing, never use a syringe inferior to 10 CC.
Risk of burst of the catheter by overpressure.

- a) Close the surgical cavity.
- b) Fix the catheter using the sterile fixing device provided.
Do not use suture for fixing the catheter.

Risk of catheter pinch-off.

- c) Connect the catheter to the perfusion line.

Always take care of the asepsis of the injection site.
Do not screw excessively the Luer connections. Fissure risk.

7 - ELIMINATION OF DEVICE AND ACCESSORIES

At the end of treatment, the catheter will be softly removed by the authorized staff.

The physician will make sure that the catheter has been entirely removed (distal end marked). The discard of the device and its accessories must follow the protocol of the discard of medical wastes DASRI (Infections Risks Care Activities Waste) related to the centers, and must not endanger the safety or health of the patient the users or any other person until complete destruction. Infectious risk by accident of exposure to blood caused by puncture or cutting accessory.

In case of inefficiency of pain treatment with our Profilis PaincAth catheter and the analgesic, a substitution treatment chosen by the medical profession will have to be performed on the patient. The catheter will be removed

Any serious incident involving the device will be notified to the manufacturer and the competent authority of the member state in which the user and / or patient is established;

Kit unitaire / Unitary set		Kit par 10 / Box of 10 sets	
P500-02.5 S/R/P		P500-02.5 QS/QR/QP	
P500-05 S/R/P	P1000-05 S/R	P500-05 QS/QR/QP	P1000-05 QS/QR
P500-07.5 S/R/P		P500-07.5 QS/QR/QP	
P500-10 S/R	P1000-10 S/R	P500-10 QS/QR	P1000-10 QS/QR
P500-15 S/R	P1000-15 S/R	P500-15 QS/QR	P1000-15 QS/QR
P500-20 S/R	P1000-20 S/R	P500-20 QS/QR	P1000-20 QS/QR
P500-25 S/R	P1000-25 S/R	P500-25 QS/QR	P1000-25 QS/QR
P500-30 S/R	P1000-30 S/R	P500-30 QS/QR	P1000-30 QS/QR

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- Place des infiltrations en peropérateurs. Injection unique vs cathéters vs systèmes à libérations prolongées. Marc Beaussier, Mourad Aissou. Département d'Anesthésie-Réanimation chirurgicale, Hôpital St-Antoine. Groupe Hospitalier Est-Parisien. APHP. Paris - Sorbonne Universités. UPMC Univ Paris 06. France. Le Congrès. Médecins. Conférence d'Actualisation. © 2014 SFAR.
est une marque déposée de / is a trademark of