disticlass médical SA

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Previously irradiated areas.	Areas near uncontrolled chronic skin lesions.
Homolateral breast cancer.	Skin metastases.
Mediastinal tumors.	Disruption of blood clotting.
A history of axillo-subclavian phlebitis.	A history of venous thrombosis.
Severe obstructive pulmonary diseases.	Patient size incompatible with device size.

10-ADVERSE SIDE EFFECTS

Percutaneous technique: Hematoma at the puncture site
Surgical technique: Hematoma in the subcutaneous lodge
 (See §15-5 OPERATING CONSEQUENCES)

11-RESIDUAL RISKS

Infectious risk.	Risk of extravasation of the drug.
Obstruction of the catheter and/or of the port.	Pneumothorax / Hemothorax / Hydrothorax.
Gas embolism.	Guide wire fragment embolism.
Superior vena cava stenosis.	Rupture and migration of the catheter.
Vein thrombosis.	Hemorrhage.
Skin ulceration and necrosis with regard to the implantable port.	Moving or turning the port.
Pericardial effusion.	Allergy to materials.
Catheter- or port-related sepsis.	Disunion of the edges of the wound.
Alteration of vessels.	Inflammation.

12-SET COMPOSITION

The DistricAth® set is delivered sterile for single use with its accessories.

	Do not re-use.		The reuse exposes patient to infection risk.
--	----------------	--	--

Implanted parts :

- 1 Implantable port DistricAth®. Unique laser-marked SN number.
- 1 Silicone Catheter or in polyurethane suited to the port exit.
- 1 Connection ring suited to the catheter diameter.

Accessories that can be integrated into the set:

• 1 Puncture needle 18 Gauges - L = 70 mm or 1 short puncture needle 21 Gauges - L = 40 mm • 1 Huber straight needle 22 Gauges - L = 25 mm • 1 90° curved Huber needle 19, 20 or 22 Gauges L = 20 to 25 mm • 1 Metallic guidewire Ø 0.018" to 0.035" L = 50 to 80 cm • 1 Tunneler Ø 1.5 to 3 mm - L = 210 mm • 1 Peel-away Introducer standard: Size = 5 to 10 French ; L = 203 mm or • 1 Peel-away Introducer short: Size = 5 to 7 French ; L = 123 mm or • 1 Valved peel-away Introducer standard : Size = 7, 9 and 10 French ; L = 196 mm • 1 Vein pick L = 50 mm Tip at 90° of 12 mm • 1 scalpel L = 143 mm - l = 11 mm • 1 Syringe 10 ml Luer L = 109 mm or • 1 Syringe 10 ml Luer Lock L = 109 mm • 1 Catheter Rinse - L = 34,5 mm	
---	--

Sterile accessories under their own EC markings with their instructions for use than can be supplied on request.

• 19, 20 gauges L=20/25mm • 1 Introducer needle 12 Gauges ; L = 80 mm • 1 cable for ECG connection L = 1 meter	
--	--

13-STERILIZATION

The DistricAth® set is delivered sterile in double packaging

	Ethylene oxide sterilized.		Do not reprocessed.
--	----------------------------	--	---------------------

	Re-sterilization can compromise the integrity of the packaging, the risk of cross-infection, device performance failure, degradation of the material and lead to toxic reactions due to residues, or even death of the patient.
--	---

14-OPENING THE KIT

Before opening this device, the nurse must have read this instruction for use in its entirety.

Carefully check the reference and expiry date of the device written on the labels of the pouch. Then check the integrity of the pouch before opening it manually to recover the blister. Then proceed to check the integrity of the blister and then open it manually in a sterile area to recover the kit elements. Take into account the visual indications when opening the packaging in order to preserve the asepsis of the device.

	Expiration date : 4 Years and 11 months Do not use this device after the expiration date written on the package.
--	---

	Do not use if the packaging is damaged. The contents of the package are sterile, as long as the package is neither visually open nor damaged, or does not show any weld breaks in the packages.
--	---

	As the devices are sterile, they must not come into contact with any non-sterile surface. Contamination of the device can cause injury, illness or even an infectious risk.
--	---

15-IMPLANTATION INSTRUCTIONS & WARNINGS

1) PRECAUTIONS AND BAD REASONABLY PREDICTABLE IMPLANTATIONS

Only qualified doctors can set up the DistricAth® implantable ports. These physicians must have demonstrated competence in the set up and use of implantable ports. They have a medical degree with specialization (Surgeon, Anesthesiologist, Radiologist, Oncologist...) Before using this device, the doctor must have read this instruction for use in its entirety.

- The implantation of implantable catheter port requires surgical asepsis and suitable equipment:

	Infectious risk for the patient.
--	----------------------------------

- Improper tightening of the puncture needle on the syringe:

	Risk of cracking of the needle by bending of the needle-syringe assembly, in particular during a difficult puncture of the subclavian vein.
--	---

- Wrong orientation of the bevel of the puncture needle:

	Risk of deformation of the guidewire spring (for example, in the case of a subclavian vein with a pronounced curvature). A blue mark on the base allows the orientation of the bevel.
--	---

- Incorrect removal of the guide wire protection cap:

	Use the guide wire pusher to remove the f from the cap so as not to damage the end of the guide wire when removing the cap.
--	---

- Blocking of the guidewire in the puncture needle :

	When catheterizing a vein (Seldinger technique), never leave the guidewire immobilized in the puncture needle in order to avoid any risk of blockage of the latter by fibrin formation.
--	---

- Dilator not properly locked in the introducer sheath.

	Ensure that the dilator base is properly locked preventing the dilator from slipping through the sheath and damaging the accessory when introducing it at skin level.
--	---

- Forced insertion of the introducer:

	For some patients, to facilitate the penetration of the introducer, provide a slight skin incision near the guidewire to avoid the risk of kinking the introducer sheath.
--	---

- Weakening of the catheter by assembly using pliers forceps:

	Catheterization with risk of rupture over time and migration into the vascular system.
--	--

- Incorrect positioning of the connection ring:

	Risk over time of catheter disconnection with migration of the catheter into the vascular network. Kit with preconnected catheter: see our technical brochure Doc. districAth®
--	---

- Rupture of the catheter by costo-clavicular pinching:

	Never implant a subclavian catheter when there is a risk of costo-clavicular pinching. The risk is rupture of the catheter and its migration into the vascular network.
--	---

- Incorrect positioning of the catheter end:

	Risk of device malfunction (no venous return by thrombosis of the catheter).
--	--

- Incorrect positioning of the port:

	Shallow lodge. Risk of necrosis of the patient's skin. Lodge too deep. Risk of incorrect location by palpation. Overstayed lodge. Risk of overturning the port.
--	---

2) IMPLANTATION SITE

Make sure you have all the necessary equipment for the implantation. Before implantation, purge the port and the catheter. The puncture in the septum using the Huber needle must be done perpendicular to the port.

The choice of the site where the implantable port will be placed may vary depending on the patient. You must select a location where the port rests on a resistant tissue (large muscle, bone surface).

The site must be adapted to the size of the port. Do not put the port too superficially to avoid any risk of necrosis of the skin. The site must be large enough to accommodate the entire port so that the scar is outside the puncture area.

The port has holes for the possible passage of stit ches to secure it on the fascia.

3) CATHETER IMPLANTATION

The catheter can be placed either by surgical technique or by percutaneous technique. Depending on the technique, the catheter can take different access ways.

When using a percutaneous introducer, it must be adapted to the diameter of the catheter.

	Two options are recommended for catheter placement:
	• Leave the catheter connected to a syringe filled with physiological saline, • Purge the catheter and clamp it at the distal end before cutting and inserting it into the vein, in order to avoid any blood rise (reduction of the risk of thrombosis).

In Europe, intravenous implantation accounts for 95% of implantable port cases.

To ensure the correct positioning of the distal end of the catheter, the latter must be at the border of the superior vena cava, atrium, controlled by an image intensifier.

The tunneling of the catheter to the port must be done manually without forcing the catheter so as not to damage it and have to start the installation again. Once, the catheter is in place in the vein, provide an additional length of catheter to make the connection with the port.

The length of the catheter should be as short as possible in order to avoid the risk of kinking of the catheter at the exit of the port. During placement, cardiac arrhythmia or tamponade related to the position of the implant may occur. This phenomenon disappears after installation.

4) CONNECTION PORT/CATHETER (See description of the connection in page 21).

Once the connection has been made and before starting the lodge, check the venous return by puncture the port with the Huber needle. Then make a lock of physiological serum or heparinized serum according to the medical protocol.

5) OPERATING CONDITIONS

The follow-up operations are generally simple.

It must be ensured that the following local complications are not present:

Hematoma next to the port or the needle site of the vein, it will subside in a few days: the room should not be used during this period if possible
Infection or inflammatory complications (redness, edema, pain...). The implantable catheter chamber should not be used until the skin is healthy again.
Skin sutures will be removed between the 7th and 10th day.

	If water-absorbable threads are used, the scar should not be wetted for 7-10 days. Risk of wound opening.
---	---

The implantable system can be used immediately after placement if the skin sutures are not directly located opposite the port.

A front chest X-ray should be prescribed. It will allow you to check the position of the catheter, the end of which should go down to the junction of the superior vena cava and the right atrium.

Fill in the card holder and the patient monitoring book.

16-USAGE & WARNINGS

1) PRECAUTIONS AND REASONABLY PREDICTABLE MISUSE

Use and maintenance may only be carried out by persons duly authorized and having received specific training.

Before using this device, the practitioner must have read this user manual in its entirety.

• The use of implantable catheter chambers requires site asepsis on the patient and unsuitable equipment:

	Infectious risk for the patient
---	---------------------------------

• Use of an unsuitable needle (coring for the septum):

	It is mandatory to use straight or curved Huber needles. The use of conventional needles would degrade the septum by coring. Risk of catheter thrombosis due to perforation.
---	--

• Incorrect needle puncture outside the septum:

	Risk of perforation of the catheter with thrombosis of the catheter. Impossibility of venous return.
---	--

• Use of an unsuitable syringe less than 10cc:

	The use of syringes of volume less than 10cc could generate overpressure a risk of bursting or rupture of the catheter with migration into the vascular system.
---	---

• Incorrect catheter rinsing:

	A rinsing with physiological saline solution shall be systematically realized at the end of each session and between two following preparations. Risk of blockage of the catheter requiring anticoagulant therapy to make the system permeable or removal of the system if treatment fails.
---	---

• MRI examination with metal accessory:

	Do not use incompatible metal accessories for NMR imaging. Risk of movement of the accessory or overheating on the patient.
---	---

2) INJECTIONS AND SAMPLES

The use of a protocol takes into account the following different steps:

• Injection of medicinal products.

• System maintenance.

• Prevention of complications related to the use of these systems (infection, rupture of the catheter, thrombosis).

These recommendations are not set in stone. Depending on the publications, adaptations could be made in particular in the area of maintenance of the implanted device (heparinized serum or not and frequency of rinsing).

Asepsis the site before any surgical gesture, skin preparation of the patient is necessary. It is carried out with an antiseptic soap from the same family as the antiseptic used in the operating room.

Access to the port must be in very specific aseptic conditions:

- The operator must wear surgical clothing, a cap, a bib and sterile gloves.
- A team of practitioners trained in the use of implantable sites.

Any injection or withdrawal is done following steps below:

- Location of the port: locate the septum of the room by palpation,
- Suturing of the Huber needle perpendicular to the skin. Crossing the septum offers resistance.

The needle should be pushed in until it stops at the bottom of the port. The correct positioning of the needle is checked by blood reflux immediately followed by rinsing by pulsed injection of a minimum of ten ml of physiological saline.

	The length of the needle must be adapted to the depth of the implantation of the chamber, and its diameter with the type of treatment. For blood samples, the first ml of blood collected must be discarded.
	Withdraw the needle while continuing to inject with the syringe to maintain positive pressure. Reduced risk of catheter thrombosis.
	When withdrawing the needle inserted into the septum, prevent a rebound effect. Risk of accidental stinging.

The use of a three-way stopcock reduces the risk of air embolism when changing syringes.

Before any injection, infusion or blood sample, the system must be rinsed with a physiological saline solution (greater than 10 ml).

Rinsing at the end of use makes it possible to make a lock of physiological serum or heparinized serum according to the medical protocol.

3) INJECTIONS

For continuous infusions, the septum puncture is usually performed with a curved Huber needle, connected to an extension tube and its 3-way stopcock. The fixation is ensured by an adhesive bandage. The occlusive dressing should allow access to the three-way stopcock. In the case of continuous infusions over several days, the needles are generally changed on average every 48 hours. During these infusions it is important to respect the established protocol.

4) SYSTEM MAINTENANCE

The purpose of maintenance is to prevent the risk of system obstruction.

The flushing rate of the system is variable according to the protocols of the establishments. Generally, the system is flushed with a physiological saline solution at a frequency of six to eight weeks in the absence of treatment.

a - Binsing volume (ISO Standard 10555-6 §4.5.4) physiological saline lock heparinized serum lock (See table p.23)

It will be essential to carry out a pulsed injection of at least 10 ml of physiological saline, followed by a lock of 5 ml of physiological saline or heparinized serum in slow injection.

NB: Doctors JL THEARD and S. ROBARD from the René Gauducheau Center in NANTES*, showed that rinsing with non-heparinized physiological saline does not increase the frequency of complications of injection site or catheter thrombosis. However, it eliminates the risk of overdose, product error and heparin intolerance. Finally, it reduces the cost and care of home maintenance by nurses, as well as the risk of infection by manipulations.

* JL THEARD, S. ROBARD, "Evolution of the protocols for using implantable chambers", Ann Fr Anesth Réanim, 1995; 14: 535

Between these infusion sessions, the patient continues his activities without constraints, except those advised against by his attending physician. Possible discomforts can be those caused by the wearing of seat belts and improper exercise.

b - Special cases

- In any case, never inject medicinal products (cytotoxic) without having checked the permeability of the system.
- If an infection is impossible, proceed to a radio control to find out the cause:
- If the catheter is thrombosed: use of anticoagulants (technique to be used only by qualified doctors).
- If the catheter is cut: removal and replacement by qualified doctors.

17-STORAGE CONDITIONS

The DistricAth implantable catheter port in its original packaging does not require any special storage conditions other than those indicated in the Symbols table in page 22.

 DistricAth implantable ports must under no circumstances be directly stored under their double sterile barrier system and must be introduced into the operating room outside the sterile area in their original packaging with the instructions for use so that the nurse and the doctor have access to the instructions for opening and installing the kit.

18-IMPLANTATION DURATION AND WITHDRAWAL

The expected lifetime from the date of implantation recommended by the manufacturer and shown on the carrier card is 5 years. It is the physician's responsibility to leave the device implanted for a longer period of time.

The duration of implantation of the DistricAth® implantable catheter port depends on the duration of treatment.

It is however recommended to remove the device according to several criteria:

- After the end of treatment depending on the psychological state of the patient.
- In the event of foreseeable non-treatment of more than 6 months following the establishment's protocol.

In addition, removal of an implantable port may be necessary, although the implant is not involved, due to side effects, such as local infection, tissue necrosis or even venous thrombosis.

The removal of this medical device must be done according to the rules of surgical asepsis with the appropriate equipment.

	When removing the implantable port, remember to remove all the elements. Make a secure withdrawal of the catheter to avoid rupture of the catheter. Risk of migration.
	The disposal of the device and of the waste must be done according to the protocol for the elimination of medical waste DASRI (Waste from Care Activities at Infectious Risks) linked to care establishments, and must not endanger safety or health of the patient, the users, or any other person until its complete destruction.
	Any serious incident occurring in connection with the device will be the subject of a notification to the manufacturer and to the competent authority of the Member State in which the user and / or patient is established. Identification of the serial number engraved on the explanted port is essential to establish the declaration. It must correspond to the serial number written on the patient monitoring book and the carrier card so that any previous declarations on the device still implanted are correct.

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est une marque déposée de / is a trademark of / es una marca depositada por / é uma marca registrada de / un marchio depositato di / is een handelsmerk van / ist ein eingetragenes Warenzeichen der